

IndoseCT

Version 20.b

Software for Calculating and Managing Radiation Dose of Computed Tomography for an Individual Patient

Manual-English

Developed by:

Choirul Anam (Diponegoro University) Fahmi R. Mahdani (Diponegoro University) Freddy Haryanto (Bandung Institute of Technology) Rena Widita (Bandung Institute of Technology) Idam Arif (Bandung Institute of Technology) Khoerun N. Syaja'ah (Bandung Institute of Technology) Geoff Dougherty (California State University Channel Islands)

Indonesia, 2020

DISCLAIMER:

This publication is based on sources and information believed to be reliable, but the authors disclaim any warranty or liability based on or relating to the contents of this publication. The IndoseCT is intended for research purposes, not for clinical applications.

COV	/ER		1					
DISC		R	2					
TAB	LE OF	CONTENTS	3					
ι.	INTR	ODUCTION	5					
	1.1	About IndoseCT	5					
	1.2	Features	5					
	1.3	Developers	8					
	1.4	List of publications	8					
	1.5	Information	10					
П.	INST	ALLATION of IndoseCT	11					
	2.1	Installation process of IndoseCT	11					
	2.2	Starting / running the IndoseCT	15					
	2.3	Creating a shortcut of IndoseCT	16					
	2.4	Uninstall of IndoseCT	19					
III.	MAII	N PART of IndoseCT	22					
IV.	PATI	ENT IMAGE	26					
	4.1	Open image	26					
	4.2	Show image	30					
	4.3	Zoom in, zoom out and move image	32					
	4.4	Image pixel value	34					
	4.5	DICOM info	35					
	4.6	Patient data	35					
v.	DETE	DETERMINATION OF THE CTDIvol						
	5.1	Input manually	38					
	5.2	Calculation	38					
		a. One slice	38					
		b. Z-axis option	42					
	5.3	From DICOM info	46					
		a. One slice	46					
		b. Z-axis option	49					
VI.	CALC	ULATION OF EFFECTIVE DIAMETER (D _{eff})	52					
	6.1	Manual calculation	53					
		a. Input D _{eff} manually	54					
		b. D _{eff} calculation from AP	54					
		c. D _{eff} calculation from LAT	55					
		d. D _{eff} calculation from LAT+AP	56					
		e. D _{eff} calculation from age	56					
	6.2	D _{eff} calculation from image	57					

		a.	Manual calculation	58
		b.	Automated calculation	59
		с.	Auto (Z-axis) calculation	64
VII.	CALCU	JLATI	ION OF THE WATER-EQUIVALENT DIAMETER (D _w)	68
	7.1	Inpu	ut the D _w value manually	68
	7.2	D _w (calculation from image	69
		a.	Manual D _w calculation	70
			D _w calculation using ellips tool	70
			D _w calculation using polygon tool	71
		b.	D _w calculation automatically (Auto)	73
		c.	D _w calculation with Auto (Z-axis)	79
VIII.	CALCU	JLATI	ION OF THE SIZE-SPECIFIC DOSE ESTIMATE (SSDE)	82
	8.1	One	e slice	82
	8.2	Z-ax	xis option	85
IX.	CALCU	JLATI	ION OF THE EFFECTIVE DOSE	87
х.	SAVE	PATII	ENT DATA	88
XI.	CALCU	JLATI	ION OF ORGAN DOSE	91
	11.1	Mo	nte Carlo data	91
	11.2	Dire	ect calculation	93
XII.	DATA	ANA	LYSIS	99
XIII.	SUMN	/IARY	,	109
XIV.	REFER	ENCE	ES	110

I. INTRODUCTION

1.1. About IndoseCT

IndoseCT 15.a is the first version of **IndoseCT**, while **IndoseCT 20.b** is the second version. **IndoseCT** was first developed in 2015. **IndoseCT** has been recorded in the Letter of Registration of Creation, Ministry of Law and Human Rights (Kementerian Hukum dan HAM), Republic of Indonesia, Number 000217029, on August 12, 2020.

IndoseCT is software to calculate and manage the radiation dose of a patient undergoing computed tomography (CT) examination. It does not only calculate the radiation output dose of CT machines (in terms of volume CT dose index, $CTDI_{vol}$), but also the individual dose received by each patient (in terms of size-specific dose estimate, SSDE) either based on the effective diameter (D_{eff}) or water-equivalent diameter (D_w). This software can also calculate radiation output and patient dose for CT equipped with tube current modulation (TCM) technique.

Until 2011, CT doses had been expressed in terms of CTDI_{vol}, whereas that metric is only able to quantify CT radiation output. The CTDI_{vol} value is influenced by almost all input parameters, such as tube voltage (kVp), tube current (mA), rotation time (s), pitch, collimation width, and others. It was noted, since its introduction, that CTDI_{vol} was not meant to indicate the radiation dose received by the patient. This is because the dose within each scanner is influenced by not only the radiation output from the CT machine, but also by the characteristics of the patient undergoing the examination, such as the size of the patient and the composition of the body part being examined. To quantify a patient's radiation dose, the size-specific dose estimate (SSDE) was introduced.

IndoseCT software estimates the radiation dose directly for each individual patient using the patient's image or by entering patient data manually. **IndoseCT** can also be used for the estimation of organ dose and effective dose. The effective dose is the dose commonly used to estimate the risk of developing cancer in the future from CT examinations. **IndoseCT** is also equipped with tools to store patient dosimetry data in a database. From this database, the software can process patient dosimetry data and easily display its graphs so that it becomes useful information for dose optimization for the use of ionizing radiation for stakeholders, whether for medical physicists, radiology doctors, hospital management or regulatory agencies.

1.2. Features

IndoseCT 15.a software features are as follows:

- Calculates the radiation output of CT devices as CTDI_{vol} for several manufacturers of CT and scanners (for various parameters of tube current (mA), tube voltage (kV), pitch, and beam collimation) (Data adopted from ImPACT 1.04)
- Calculates the radiation output of CT devices as CTDI_{vol} for scanners equipped with tube current modulation (TCM) or automatic exposure control (AEC) techniques (Anam et al. Int J Rad Res. 2018; 16(3): 289-297).
- 3. Calculates the effective diameter (D_{eff}) for variations in age, lateral (LAT) diameter, anterior-posterior (AP) diameter, and combined LAT+AP). This feature is similar to previous available calculators. The user only inputs patient data (in the form of patient

age, and LAT diameter AP diameter, or a combination of LAT+AP) and the D_{eff} value is immediately obtained (Data adopted from AAPM TG 204)

- 4. Calculates the effective diameter (D_{eff}) from a CT image for 2D or 3D, both manually and fully automatically. Manual calculation means that the user performs measurements manually with the **line** tool on the patient image. Fully automatic calculation means that the D_{eff} value of the patient image is determined without intervention from the user (Anam et al. Adv Sci Eng Med. 2015; 7: 892-896). In fully automatic calculations, the user can also select the diameter at the maximum position, the diameter at the center position, or directly use the cross-sectional area of the patient image (Anam et al. Atom Indonesia. 2017; 43(1): 55-60).
- 5. Calculates the water-equivalent diameter (D_w) of CT images for 2D or 3D (manually and automatically) (Anam et al. J Appl Clin Med Phys. 2016; 17(4): 320-333)
- Calculates the D_w value for a truncated image that is cut off at the edges. Truncated images are common in clinical applications. IndoseCT is equipped with a correction factor for truncated images (Anam et al. Radiat Prot Dosim. 2017; 175(3): 313-320).
- Calculates and displays the D_{eff} and D_w profiles along the longitudinal axis for all images. D_{eff} and D_w calculations for all 3D CT images may require a relatively long computational time. **IndoseCT** also allows the user to select the number of images to be calculated or use certain intervals (eg images 1, 10, 30, and so on).
- 8. Calculates patient dose in terms of size-specific dose estimate (SSDE) from the $CTDI_{vol}$ and D_{eff} (or D_w) values entered by the user, or directly from CT images for 2D or 3D (manually and automatically) (Anam et al. .J Phys Conf Ser. 2016;694:012030).
- 9. Calculates the total radiation dose as the metric of dose-length product (DLP), either as standard or corrected by patient size.
- 10. Calculates organ doses based on patient size taken directly from patient CT images for various protocols (data and equations adopted from Sahbaee et al. Med Phys. 2014;41(7):072104).
- 11. Calculates the effective patient dose based on patient size for various examination protocols (data and equations adopted from Sahbaee et al. Med Phys. 2014; 41(7): 072104).
- 12. Extracts some patient data from the DICOM info and stores it in the database. This data storage and processing will be very beneficial for the institution, so that it can manage doses and other data locally, and can take strategic steps related to the application of CT in the institution. This data is relatively small and is stored separately from patient images which require very large storage media.
- Analyzes patient dosimetry data and displays it as needed. There are many data options that can be displayed, such as CTDI_{vol}, DLP, D_{eff}, D_w, SSDE profiles, as well as relationships between metrics such as the relationship between D_{eff} and SSDE, or the relationship between CTDI_{vol} and SSDE, and so on.

The development features of **IndoseCT 20.b** are as follow:

- 1. **IndoseCT 20.b** can be run without having to install the main program (Matlab). Users can directly install or uninstall **IndoseCT** from the computer easily.
- 2. The display of IndoseCT 20.b is very user-friendly so that it is easier to use.
- 3. **IndoseCT 20.b** is faster at performing calculations than the previous version.

- 4. Access options to open files using folders, while earlier versions only had a file-based option. Opening an image within a folder is faster than as a file. In this case there are no other files in the folder other than the patient image file.
- 5. Sample images have been added in the form of anthropomorphic phantom images from the base of the pelvis to the top of the head. These samples are very useful for practice in using **IndoseCT**.
- 6. It has been equipped with several window options to display images with higher contrast between objects, such as soft tissue, bone, and so on.
- 7. To move from the image of one slice to another, the arrow keys on the keyboard can be used in this version so that it is more practical.
- 8. The image can be zoomed-in, zoomed-out, and shifted using the keyboard.
- 9. DICOM info is easy to access. Users sometimes want to access such information, for example information about the field of view (FOV).
- 10. A database of scanner types and CTDI_{vol} values for several manufacturers and scanners has been added (Data adopted from ImPACT 1.04 and WAZA-ARI).
- 11. There is an option to display tube current profiles, CTDI_{vol}, and SSDE (Anam et al. Information. 2017; 20(1): 377-382)
- 12. It has been equipped with an additional option to adjust the CTDI_{vol} value in the TCM technique. For example, for the image taken with the TCM technique the CTDIvol value does not fluctuate with the current, but we have added an option to adjust the CTDIvol value.
- In some cases CTDI_{vol} value, may not be contained in the DICOM info, so then other parameters (such as tube current, tube voltage, etc.) can be taken automatically from the DICOM info to calculate the CTDI_{vol} value. With this facility, CTDI_{vol} calculation becomes faster.
- 14. D_{eff} calculations have been complemented by options for correcting for the presence of lung (As introduced by Mihailidis et al. Br J Radiol. 2020; 93: 20200473) and bone.
- 15. For the calculation of D_w, a new algorithm has been developed to detect the patient's body as a whole, even though there are several separate parts, for example when there are two patient arms on the thoracic image (Anam et al. J Appl Clin Med Phys. 2021;1-11). This is different from the previous version which only uses the largest patient object selection, so that when there are two patient arms in the thoracic image, only one part is segmented, namely the thorax, while the two arms are not segmented.
- 16. For D_w calculations, there is an option to calculate D_w from the entire image (without segmentation) and there is a further option to remove the patient table from the image automatically (Anam et al. Radiat Prot Dosim. 2019; 185(1): 42-49)
- 17. For 3D (Z-axis) options in the calculation of D_{eff} and D_w , it is equipped with **Regional** options, so the user can determine the calculated slice, for example calculating the value from the 51st slice to the 100th slice. This is useful for calculating organ doses (Anam et al. J Biomed Phys Eng . 2021)
- 18. There is a new size-conversion value option (from CTDIvol to SSDE) for head examinations (Data adopted from AAPM TG 293).
- 19. A system has been equipped for calculating dose distribution within a patient (dosemap) and calculating organ doses (Anam et al. J X Ray Sci Med. 2020; 28: 695-708), although the contouring process should be still done manually by the user.

- 20. Patient data stored in the database can be deleted easily. It can also be exported to Microsoft Excel easily.
- 21. Data stored in databases such as institutional data, vendors, and scanner types, can be mixed with data from various scanners from one institution (hospital) or several institutions for comparison. **IndoseCT** will differentiate the data and process it according to the needs and preferences of users.
- 22. The **IndoseCT** can display graphs (e.g. SSDE and D_w relationships), and it is equipped with trendline options such as linear equations, quadratic equations, polynomial equations, exponential equations and others.
- 23. The displayed graph is equipped with the option to display the average value and standard deviation. The x-axis and y-axis values can be adjusted as needed. Graphics can also be shifted for a more optimal display.
- 24. The displayed graph can be saved in various image formats (such as jpeg and bitmap) or exported to Microsoft Excel.

1.3. Developers

- Choirul Anam (*Diponegoro University*)
- Fahmi R. Mahdani (*Diponegoro University*)
- Freddy Haryanto (*Bandung Institute of Technology*)
- Rena Widita (Bandung Institute of Technology)
- Idam Arif (*Bandung Institute of Technology*)
- Khoerun N. Syaja'ah (*Bandung Institute of Technology*)
- Geoff Dougherty (California State University Channel Islands)

1.4. Publications

For those who want to delve deeper into the scientific aspects of **IndoseCT**, please refer to the following list of publications:

- 1. Anam C, Haryanto F, Widita R, Arif I. *Automated estimation of patient's size from 3D image of patient for size specific dose estimates (SSDE)*. **Adv Sci Eng Med.** 2015; 7(10): 892-896.
- Anam C, Haryanto F, Widita R, Arif I, Dougherty G. The evaluation of effective diameter (D_{eff}) calculation and its impact on size-specific dose estimate (SSDE). Atom Indonesia. 2017; 43(1): 55-60.
- 3. Anam C, Haryanto C, Widita R, Arif I, Dougherty G. *Automated calculation of water-equivalent diameter (D_w) based on AAPM report TG. 220.* **J Appl Clin Med Phys.** 2016; 17(4): 320-333.
- 4. Anam C, Mahdani FR, Dewi WK, Sutanto H, Triadyaksa P, Haryanto F, Dougherty G. *An improved method for automated calculation of the water-equivalent diameter for estimating size-specific dose in CT*. J Appl Clin Med Phys. 2021; 22(9): 313–323.
- 5. Anam C, Haryanto F, Widita R, Arif I, Dougherty G. *A fully automated calculation of size-specific dose estimates (SSDE) in thoracic and head CT examinations*. **J Phys Conf Ser.** 2016; 694: 012030.
- 6. Anam C, Haryanto F, Widita R, Arif I, Dougherty G. *The size specific dose estimates (SSDE)* for truncated computed tomography images. Radiat Prot Dosim. 2017; 175(3): 313-320.

- 7. Anam C, Haryanto F, Widita R, Arif I, Dougherty G, McLean D. *The impact of patient's table on size-specific dose estimates (SSDE)*. Australas Phys Eng Sci Med. 2017; 40(1): 153-158.
- 8. Anam C, Arif I, Haryanto F, Widita R, Lestari FP, Adi K, Dougherty G. A simplified method for the water-equivalent diameter calculation to estimate patient dose in CT examinations. Radiat Prot Dosim. 2019; 185(1): 42–49.
- 9. Anam C, Haryanto F, Widita R, Arif I, Dougherty G, McLean D. Volume computed tomography dose index (CTDIvol) and size-specific dose estimate (SSDE) in tube current modulation CT. Int J Rad Res. 2018; 16(3): 289-297.
- 10. Anam C, Fujibuchi T, Toyoda T, Sato N, Haryanto F, Widita R, Arif I, Dougherty G. *A simple method for calibrating pixel values of the CT localizer radiograph for calculating water-equivalent diameter and size-specific dose estimate*. **Radiat Prot Dosim**. 2018; 179(2): 158–168.
- 11. Anam C, Haryanto F, Widita R, Arif I, Dougherty G. *Profile of CT scan output dose in axial and helical modes using convolution*. **J Phys Conf Ser.** 2016; 694: 012034.
- 12. Anam C, Haryanto F, Widita R, Arif I, Dougherty G. *The profile of patient dose along the longitudinal axis in CT using tube current modulation (TCM)*. **Information (Japan).** 2017; 20(1B): 377-382.
- 13. Anam C, Haryanto F, Widita R, Arif I, Dougherty G. *Estimation of eye radiation dose during nasopharyngeal CT examination for an individual patient*. Information (Japan). 2016; 19(9B): 3951-3962.
- 14. Anam C, Haryanto F, Widita R, Arif I, Dougherty G. A size-specific effective dose for patients undergoing CT examinations. J Phys Conf Ser. 2019; 1204: 012002.
- 15. Anam C, Adhianto D, Sutanto H, Adi K, Ali MH, Rae WID, Fujibuchi T, Dougherty G. *Comparison of central, peripheral, and weighted size-specific dose in CT.* J X-ray Sci Technol. 2020; 28: 695–708.
- 16. Anam C, Dewi WK, Masdi M, Haryanto F, Fujibuchi T, Dougherty G. *Investigation of eye lens dose estimate based on AAPM report 293 in head computed tomography*. J Biomed Phys Eng. 2021: 1-10.
- 17. Anam C, Fujibuchi T, Haryanto F, Widita R, Arif I, Dougherty G. *An evaluation of computed tomography dose index measurements using a pencil ionisation chamber and small detectors*. J Radiol Prot. 2019; 39: 112–124.
- 18. Adhianto D, Anam C, Sutanto H, Ali MH. *Effect of phantom size and tube voltage on the size-conversion factor for patient dose estimation in computed tomography Examinations*. **Iran J Med Phys**. 2020; 17: 282-288.

The following is a list of several publications on studies conducted using IndoseCT:

- 1. Anam C, Budi WS, Adi K, Sutanto H, Haryanto F, Ali MH, Fujibuchi T, Dougherty G. *Assessment of patient dose and noise level of clinical CT images: automated measurements*. **J Radiol Prot**. 2019; 39: 783–793.
- 2. Fahmi A, Anam C, Suryono, Ali MH. *The size-specific dose estimate of paediatric head CT examinations for various protocols*. **Radiat Prot Dosim**. 2020; 188(4): 522–528.
- 3. Fahmi A, Anam C, Suryono, Ali MH, Jauhari A. *Correlation between age and head diameters in the paediatric patients during CT examination of the head*. **Pol J Med Phys Eng**. 2019; 25(4): 229-235.

- 4. Utami MSN, Sutanto H, Anam C. Effect of contrast agent administration on size-specific dose estimates (SSDE) calculations based on water equivalent diameter in CT head examinations. International Journal of Scientific Research in Science and Technology. 2021; 8(3): 563-571.
- 5. Dewi WK, Anam C, Hidayanto E, Nitasari A, Dougherty G. *The effective and water-equivalent diameters as geometrical size functions for estimating CT dose in the thoracic, abdominal, and pelvic regions*. **Pol J Med Phys Eng**. 2021.
- 6. Dewi WK, Anam C, Hidayanto E, Wati AL, Dougherty G. *Correlation between anterior posterior and lateral dimensions and the effective and water-equivalent diameters in axial images from head computed tomography examinations*. **Radiat Prot Dosim**. 2021.
- 7. Nitasari A, Anam C, Budi WS, Wati AL, Syarifudin S, Dougherty G. *Comparisons of water-equivalent diameter measured on images of abdominal routine computed tomography with and without a contrast agent*. **Atom Indonesia**. 2021: 1-5.
- 8. Ali MH, Anam C, Haryanto F, Dougherty G. *The comparison of size-specific dose estimate in CT examination based on head and body PMMA phantom*. **Journal of Physics and Its Applications**. 2018; 1(1): 1-6.
- Wati AL, Anam C, Budi WS, Nitasari A, Syarifudin, Satoto B, Jannah M, Sepsatya F. Patient doses based on the acceptable quality dose on thoracic CT examination. International Journal of Latest Engineering Research and Applications. 2021; 6(6): 17-22.

1.5. Information

For further information about **IndoseCT 20.b** software, please contact email: <u>anam@fisika.fsm.undip.ac.id</u> or <u>anamfisika@gmail.com</u>

II. INSTALLATION of IndoseCT

2.1 Installation process of IndoseCT

Before installing IndoseCT, first make sure that the **IndoseCT** installation file has been obtained and downloaded to your computer. To get the **IndoseCT** installation file and password, you can contact the email given in the Information section (Chapter I). The **IndoseCT** installation file size is relatively small, only around 130 MB. There are two choices of **IndoseCT** installation files, namely for **32-bit** and **64-bit** systems. The **IndoseCT** installation file for a **64 bit** system is shown in **Figure 1**.



Figure 1. IndoseCT installation file. In this case the file name is IndoseCTSetup.v20b.21 64 bit.

To install **IndoseCT**, double click on the file. Next, information will appear that **Microsoft Defender** has detected that this file is an unrecognized application that has the potential to harm your computer (as shown in **Figure 2**). This is not a problem. Click **More info** to get a box as shown in **Figure 3**, and click **Run anyway** to continue the installation process.

If you want to continue to install, select **Yes** (from **Figure 4**), and enter the **password**. Next, click **OK** to continue

After that, **IndoseCT Setup** will appear. To continue click **Next (Figure 5)**. The **Destination Folder** option will appear as shown in **Figure 6**.



Figure 2. Information shows **Microsoft Defender** detects that this file is an unrecognized application that has the potential to harm your computer. This is not a problem. Continue by clicking **More info**.



Figure 3. Information shows that this application is **IndoseCTSetup.v20b.23 64 bit**. Click **Run anyway** to continue.

User Account Control Do you want to allow this app to make changes to your device?	×
Uninstall or change an application Verified publisher: Microsoft Windows	
Yes No .	

Figure 4. The dialog contains options whether you want to allow this software to be installed on your computer or not.



Figure 5. IndoseCT Setup dialog. Select Next button to continue the installation process.



Figure 6. Selection of the destination folder to place the **IndoseCT** file to be installed on the computer.

The **Destination Folder** by default is located in **C:\Program Files\IndoseCT**. If you agree to this **Destination Folder**, the installation process will continue. (Or if you want to place it in a specific folder, select **Browse**, enter the desired **Destination Folder**. Click **Install**, and the installation progress dialog will appear as shown in **Figure 7**. Wait a few moments for the installation process to complete. The time of the installation process takes no more than 5 minutes. When finished, a **Completed** notification will appear. Then click **Next**, and the **Completing IndoseCT Setup** dialog appears. Before choosing **Finish**, you can choose the **Run IndoseCT** option to open **IndoseCT** opened immediately, after clicking the **Finish** button.

At this point, the **IndoseCT** installation is completed. IndoseCT is ready to be used to calculate patient doses. The first screen of **IndoseCT** is shown in **Figure 8**.

This PC	MATLAB R2014	11/200	
Recycle Bin	Media Player Classic	🕼 IndoseCT Setup — 🗆 🗙	
Control Panel	Microsoft Edge Microsoft	Installing Please wait while IndoseCT is being installed.	9 a
Adobe Reader XI	paintmet New shortcut	Show details	
Firefox	Polyspace Bug Find		
imquest	Polyspace Visual Studio Code Prov Visual Studio	Nullsoft Install System v3.06.1	
ari Jitsi	PyCharm WirFi- Communi Shortcut		

Figure 7. IndoseCT installation progress dialog.

G IndoseCT v20b									-	٥	×
	\$ 2.	Windowing: Soft Tis	sue 🗸 360 60						Phantom	HEAD	~
0 -					ID Name Age Sex CTDJ Meth Calc	- Unspecified Vol Diamete	r SSDE Or	Exam Date Institution Scanner Protocol	16/09/2021 Manufacturer-Mod	el V	
200 - 300 -					One Para Ma Scca Vol	slice ameter nufacturer nner :age (kV)	Siemens Emotion 110	Output mAs Effective CTDI _w (n	mAs	0 0 0 0	
400 -					Rot Pite Col	e Current (mA) ation Time (s) h imation (mm)		DLP (m	iy-cm)	0	
500 - 600 - -100	0	100 200	300 400	500	600 200	Calculate	Get Info			Next]
Sort Images	< 0/0 > 0	Go to slice									

Figure 8. Main view of IndoseCT 20.b.

2.2. Starting / Running the IndoseCT

After installation, **IndoseCT** is ready to be used anytime. Left-click the **Start** button at the lower left corner, and click **IndoseCT**.

Zoom	GOM Player MATLAB R2014	in the second			1-32			
≡	OM Player	Create			Play			
	Groove Music	Sunday				5 0		
	1	8	6					-T-IN
	imquest	0	Mail		Xbox Console		Photos 📑	
	indoQCT v				Sex			
L			0		ENDS			
	Intel® Graphics Command Center		Office	41	FRIE	Solitaire		
	Intel® Graphics Control Panel	Explore			Lenovo			
	J				1			
	🦰 JetBrains 🗸 🗸		•					
	🤞 Jitsi	Microsoft Store	_	Microsoft Edge	Vantage			
8	К	Fair	Play music at					
n	K-Lite Codec Pack ~	25° ^{33°} 25°	home or on the go.	S				
	KMSAuto Net 2015 1.4.0 ~	Jakarta	Spotify 🕒	Skype	6			
2	L							
0	💼 Lenovo Hotkeys		۲		IndoseCI			
()	Lenovo Vantage	News 🛄	Paint 3D					
	in LinkedIn							
		C) 🔒 🖥		x w	I 🔁		🤙 26°C Kabut ∧ 🗈 📴 08/08/2021 💀

Figure 9. Opening IndoseCT via the Start menu button.

2.3. Creating a shortcut for IndoseCT

To be able to open **IndoseCT** quickly, an **IndoseCT** shortcut can be created and displayed on the desktop. To create an **IndoseCT** shortcut, right-click on the desktop, select **New**, then select **Shortcut** as shown in **Figure 10**.



Figure 10. The display showing how to create a shortcut forf IndoseCT.

Next, select **Browse** for the folder where **IndoseCT** was placed during installation. If during installation you left the default option (**C:\Program Files\IndoseCT**), then **IndoseCT** is in **C** and then in **Program Files (Figure 11)**.

This PC	MATLAB R2014 Media Player Classic	GOM Player		17		*		
Control Panel Adobe Reader XI	Microsoft Edge paintanet	Microsoft Teams	 x. Create Shortcut What item woul This wizard helps you Internet addresses. Type the location of 	Browse for Files or Folders Select the target of the shortcut below:	×	s, computers, or Browse		
Firefox	Polyspace Bug Find	Visual Studio Code	Click Next to contine	Common rues imagea indeseCT intel Maix New Folder OK	Cancel	Next Cancel		
a Jitsi	PyCharm Communi	Zoom						

Figure 11. IndoseCT folder is located in C:\Program Files

Look for the **IndoseCT.exe** file (as shown in **Figure 12**), select it and click **OK**. Then the exe-file is ready to be the **IndoseCT** shortcut. Select **Next** to continue, and then click **Finish** as shown in **Figure 13**.

This PC	MATLAB R2014		
Recycle Bir	Media Player Classic	K K	
Control Panel Adobe Reader XI	Microsoft Edge Microsoft reams paintant RUNDLG RUNDLG Polyppine Visual Studio Bug Finda.	Browse for Files or Folders X Select the target of the shortcut below: x, computers, or This wizard helps yot Internet addresses. api-ms-win-crt-utility-11-1-0. Type the location of api-ms-win-crt-utility-11-1-0. Click Next to contin base_library Click Next to contin timbaseCT Window News Endow timbasecT Window News Endow timbasecT Window News Endow timbasecT	
a imquest a Jitsi	Polyspace Code Prov. PyCharm PyCharm	Net Cancel	

Figure 12. IndoseCT.exe file.

This PC	MATLAB R2014	GOM Player	APRIL CONTRACT	
Recycle Bin	Media Player Classic	IndeQCT	K K	
Control Panel	Microsoft Edge	Microsoft Teams	What item would you like to create a shortcut for? This wizard helps you to create shortcuts to local or network programs, files, folders, computers, or	
Adobe Reader XI	paint-net	RUNDLG	Type the location of the item: Type the location of the item: Browse	
Firefox	Polyspace Bug Find	Visual Studio Code	Click Next to continue.	
imquest	Polyspace Code Prov	Wi-Fi - Shortcut	Next Cancel	
E Jitsi	PyCharm Communi	Zoom		

Figure 13. The location of the **IndoseCT.exe** file has been selected and is ready to be displayed as a shortcut on the desktop.

You can sort the shortcuts on the desktop, for example by name, by right-clicking on the desktop, then **Sort by**, then selecting **Name**. At that time, all the shortcuts on our desktop have been sorted by name. An example of the **IndoseCT** shortcut on the desktop, is shown in **Figure 14**.



Figure 14. IndoseCT shortcut on computer desktop. To enter **IndoseCT**, all you have to do is to double click on the shortcut.

2.4. Uninstall the IndoseCT

The process of uninstalling IndoseCT is also very easy. First, left-click the Start button and search for IndoseCT. If it is visible, then right-click IndoseCT (Figure 15). Then select Uninstall. You can also go directly to Control Panel Home by clicking the Control Panel shortcut (Figure 16). After that, click Programs and Features (Figure 17).



Figure 15. The process of uninstalling **IndoseCT**. Click the **Start button** (in this example at the lower left corner of the desktop), then search for **IndoseCT**. If you have found it, **right-click** it, and the **Uninstall** option will appear.



Figure 16. Shortcut Control Panel.

🛃 All Control Panel Items				– 0 ×
← → · · ↑ 🖾 > Control Panel > Al	I Control Panel Items >			م ي ا
Adjust your computer's settings				View by: Large icons 🔻
Administrative Tools	AutoPlay	Backup and Restore (Windows 7)	Color Management	🗿 Credential Manager
Date and Time	Default Programs	Device Manager	Devices and Printers	Ease of Access Center
File Explorer Options	File History	Fonts	Indexing Options	Part Internet Options
	Mail (Microsoft Outlook)	🥏 Mouse	Network and Sharing Center	Phone and Modem
Power Options	Programs and Features	kecovery	Region	RemoteApp and Desktop Connections
Y Security and Maintenance	Sound	Speech Recognition	Storage Spaces	Sync Center
System	Taskbar and Navigation	Troubleshooting	🝇 User Accounts	Windows Defender Firewall
Windows Mobility Center	Work Folders			

Figure 17. To uninstall IndoseCT, double-click on Programs and Features.

The **Control Panel Home** dialog will open as shown in **Figure 18**. Left-click **IndoseCT**, then **Uninstall/Change** will appear at the top. If you want to continue the uninstall process, then click **Uninstall/Change**. Next, there will be a notification whether you want to continue this process or not (**Figure 19**). If so, select **Uninstall**.

 ← → · · · ↑ ○ · Control Panel · › All Control Panel · · › All Control Panel · · · · · · · · · · · · · · · · · · ·	eader XI (11.0.23) tings	click Uninstall, Change, or Repair. Publisher Adobe Systems Incorporated Advanced Micro Devices, Inc.	Installed On 29/08/2019	Size		~	Ğ		Q==	م
Control Panel Home Unins View installed updates To unins off Turn Windows features on or off Crganice - Name Mano Set Organice - Name Mano Set Mano Set Ma	stall or change a program stall a program, select it from the list and then Uninstall/Change eader XI (11.0.23) tings yer	click Uninstall, Change, or Repair. Publisher Adobe Systems Incorporated Advanced Micro Devices, Inc.	Installed On 29/08/2019	Size						
View installed updates To unins Turn Windows features on or off Organize Name Adobte Organize Mame Marce Imag	stall of change a program stall a program, select it from the list and then Uninstall/Change eader XI (11.0.23) tings yer	click Uninstall, Change, or Repair. Publisher Adobe Systems Incorporated Advanced Micro Devices, Inc.	Installed On 29/08/2019	Size						
View installed updates To uninit off Turn Windows features on or off 300 Play 300 M Pla	stall a program, select it from the list and then Uninstall/Change eader XI (11.0.23) tings type	click Uninstall, Change, or Repair. Publisher Adobe Systems Incorporated Advanced Micro Devices, Inc.	Installed On 29/08/2019	Size						
Turn Windows features on or off Organize Name Adobe Re JB Adobe Re JE IndoseCT LE IndoseCT LE IndoseCT LE IndoseCT LE IndoseCT LE IndoseCT LE IndoseCT LE IndoseCT LE IndoseCT LE INFO Microsoft	eader XI (11.0.23) tings yer	Publisher Adobe Systems Incorporated Advanced Micro Devices, Inc.	Installed On 29/08/2019	Size						
off Organize • Name Made Re Marker Ma	eader XI (11.0.23) tings	Publisher Adobe Systems Incorporated Advanced Micro Devices, Inc.	Installed On 29/08/2019	Size						
Name Adobe R	eader XI (11.0.23) tings yer	Publisher Adobe Systems Incorporated Advanced Micro Devices, Inc.	Installed On 29/08/2019	Size					8== 👻	
Adobe Re jj Adobe Re jj Adobe Re i GOM Pia or mqueze i mqueze i mqueze i mose A i mose	eader XI (11.0.23) tings tyer	Adobe Systems Incorporated Advanced Micro Devices, Inc.	29/08/2019		Version					
MADE RE MAN Set M	itings Iver	Advanced Micro Devices, Inc.	2 7 () W V () () 7	100 M/P	11.0.22					
GOM Pia GOM Pia Gom Pia Indusect	iyer	Advanced Micro Devices, inc.	20/02/2021	190 1910	2019 1004 2240 42895					
Industry		GOM & Company	30/03/2021		2 3 30 5289					
Initiation I		time timpersity	50/05/2021		///					
Intri Mirris Kiris Kiris Kiris Kiris Kiris Kiris Kiris Kirise Kirise Martaß Martaß Microsoft Mi	(Fahmi Rosdiansvah Mahdani	06/08/2021	308 MB	10					
itsi itsi itsi itsi itsi itsi itsi its	Annagement Engine Components	Intel Corporation	30/ 03/ 2021	107 110	1024.12.0.1140					
K-Lite Co E (MKSAuto Liencov MASAuto MartLAB MATLAB Microsofi Micro		Jitsi	24/03/2020	175 MB	2.10.5550					
KMSAuto KmSAuto Kenvo V MartLaB MartLaB Microsoft Microso	odec Pack 6.5.0 (Full)		21/08/2019	48.2 MB	6.5.0					
Lenovo V A MATLAB A MATLAB Microsoft	o Net 2015 1.4.0		01/05/2021	9.05 MB						
MATLAB MATLAB MATLAB Microsoft Micros	Vantage Service	Lenovo Group Ltd.	03/06/2021	25,9 MB	3.7.19.0					
MATLAB Microsoft	R2014b (32-bit)	The MathWorks, Inc.	30/03/2021		8.4					
Microsoft	Runtime 9.3	MathWorks	30/03/2021		9.3					
C Microsoft Microsoft Microsoft Microsoft Microsoft	ft Access 2019 - en-us	Microsoft Corporation	06/08/2021		16.0.14228.20226					
Microsoft Microsoft Microsoft Microsoft Microsoft Microsoft	ft Edge	Microsoft Corporation	07/08/2021		92.0.902.67					
 Microsoft Microsoft Microsoft Microsoft 	ft Excel 2019 - en-us	Microsoft Corporation	06/08/2021		16.0.14228.20226					
Microsoft Microsoft Microsoft	ft OneDrive	Microsoft Corporation	03/08/2021	188 MB	21.139.0711.0001					
C Microsoft	ít Outlook 2019 - en-us	Microsoft Corporation	06/08/2021		16.0.14228.20226					
1 Missen	ít PowerPoint 2019 - en-us	Microsoft Corporation	06/08/2021		16.0.14228.20226					
	ít Project - en-us	Microsoft Corporation	06/08/2021		16.0.14228.20226					
1 Microsoft	ít Publisher 2019 - en-us	Microsoft Corporation	06/08/2021		16.0.14228.20226					
🗊 Microsoft	ít Teams	Microsoft Corporation	15/07/2021	114 MB	1.4.00.16575					
Microsoft	ít Update Health Tools	Microsoft Corporation	06/08/2021	1,09 MB	2.82.0.0					
1 Microsoft	ít Visio - en-us	Microsoft Corporation	06/08/2021		16.0.14228.20226					
Microsoft	ít Visual C++ 2005 Redistributable	Microsoft Corporation	24/08/2019	4,84 MB	8.0.61001					
Microsoft	rt Visual C++ 2008 Redistributable - x86 9.0.3	Microsoft Corporation	24/08/2019	10,1 MB	9.0.30729.6161					
Hicrosoft	ft Visual C++ 2013 Redistributable (x64) - 12	Microsoft Corporation	30/03/2021	20,5 MB	12.0.30501.0					
		er en la			****					
		IB								

Figure 18. Control Panel Home dialog for uninstalling IndoseCT.

Programs and Features									a.	_	٥	×
Control i	aner > All Control Panel items > Program	s and realures						~	0			þ
Control Panel Home	Uninstall or change a progra	m										
View installed updates	To uninstall a program, select it from t	he list and then	click Uninstall Change or Ren	ur.								
Turn Windows features on or off			ener ennsten, energe, er rep									•
	Organize • Uninstall/Change										8== *	Ø
	Name	🛞 IndoseCT L	Ininstall		- 0	×						
	🖊 Adobe Reader XI (11.0.23)											
	AMD Settings	Doministrali Ind	IOSEC I			Vanis -	2349.42886					
	🔁 GOM Player	Kenove ind	used i from your computer.				9					
	imquest											
	IndoseCT	IndoseCT w	I be uninstalled from the following	folder. Click Uninstall to s	tart the							
	Management Engine Componer						1140					
	KMSAuto Not 2015 1 4 0	Uninstalling	from: C:\Program Files\Indose	21								
	Lenovo Vantage Service											
	MATLAB R2014b (32-bit)											
	MATLAB Runtime 9.3											
	Microsoft Access 2019 - en-us						.20226					
	C Microsoft Edge						7					
	Microsoft Excel 2019 - en-us						.20226					
	 Microsoft OneDrive 						1.0001					
	🚺 Microsoft Outlook 2019 - en-us	Nullsoft Install :	ivstem v3.06.1				.20226					
	Microsoft PowerPoint 2019 - en-us			11-1	-	Para and	.20226					
	Microsoft Project - en-us			Unins		Jancel	.20226					
	Microsoft Publisher 2019 - en-us		Microsoft Corporation	06/ 06/		10.0.1422	rd.20226					
	Microsoft Teams		Microsoft Corporation	15/07/2021	114 MB	1.4.00.16	575					
	Microsoft Update Health Tools		Microsoft Corporation	06/08/2021	1,09 MB	2.82.0.0						
	Microsoft Visio - en-us		Microsoft Corporation	06/08/2021	4.04.140	16.0.1422	8.20226					
	Microsoft Visual C++ 2008 Pedistributa	ble - v86 9 0 ?	Microsoft Corporation	24/06/2019	4,04 MB	0.0.01001	6161					
	Microsoft Visual C++ 2008 Redistributa	ble (x64) - 12	Microsoft Corporation	30/03/2021	20.5 MB	12.0.3050	1.0					
	July and the second second			30/03/2021	20,5 110	12.0.5050						
	Fahmi Rosdiansyah Mahd	ani Size: 308 N	1B									

Figure 19. Dialog for notification whether you continue the uninstall process or not.

A dialog on the progress of the uninstall process appears, and when it is finished, a **Completed notification** will appear as shown in **Figure 20**. Then select **Close**. At this point, the process of uninstalling **IndoseCT** is complete.

This PC	MATLAB R2014	Wi-Fi - Shortcut	-4-	and the		A KIN
Recycle Bin	Media Player Classic	Zoom	IndoseCT Uninstell	×		NAME OF A
Control Panel	Microsoft Edge		Uninstallation Complete Uninstall was completed successfully. Completed	6	ý -	
Adobe Reader XI	paint.net		Show details			
Firefox	Polyspace Bug Find					
imquest	Polyspace Code Prov		Nulsoft Install System v3.06.1	Close Cancel		
IIIII	PyCharm Communi, Visual Studio Code					

Figure 20. Notification that the uninstall process is complete.

III. MAIN PART of IndoseCT

In the previous chapter, we discussed how to install and run **IndoseCT**. The main view of **IndoseCT 20.b** after opening was shown in **Figure 8**. This chapter will discuss the parts of **IndoseCT** in general.

The <u>first part</u> of **IndoseCT** is the main part to calculate patient size, radiation dose and analyze the data that has been obtained. This section consists of 5 tabs, namely **CTDIvol**, **Diameter**, **SSDE**, **Organ**, and **Analyze** (Figure 21).

The **CTDI**_{vol} **tab** is a section for calculating , obtaining or entering the CTDI_{vol} value. CTDI_{vol} itself is a quantity that describes the output dose of the CT machine.



Figure 21. The main view of **IndoseCT**, which consists of several tabs to calculate patient size, radiation dose and analyze data.

The **Diameter tab** is a section for calculating, obtaining or simply entering a patient's diameter value. In general, two diameters are used, namely the effective diameter (D_{eff}) and the water-equivalent diameter (D_w).

The **SSDE tab** is a section for calculating size-specific dose estimate (SSDE) values. SSDE is a quantity that describes the dose received by the patient. SSDE is measured from the device output dose (CTDI_{vol}) and patient size characteristics (diameter). Therefore, to calculate the SSDE value, the previous two tabs must be completed (**CTDI_{vol}** and **Diameter tabs**)

The **Organ tab** is a section for calculating radiation doses to an organ. Calculation of the dose of organs is based on the SSDEs. In order to obtain an organ dose value, the previous three tabs have to be completed (**CTDI**_{vol}, **Diameter**, and **SSDE tabs**). However, in **IndoseCT** organ doses are not stored in the database, only CTDI_{vol}, diameter and SSDE values are stored.

The Analyze tab is a section for analyzing the results that have been stored in the database. Various parameters can be connected and displayed in graphical form, for example the relationship between D_w and SSDE, and so on.

Detailed discussion of each tab will be discussed in later chapters.

The <u>second part</u> of IndoseCT is a screen for displaying patient images (Figure 22), including segmentation results and various attempts to get the diameter displayed on this screen. This screen can be enlarged or minimized by dragging the right border of the image to the right or to the left.

a IndoseCT v20b I 🖿 🖿 略 線 👰 🖥 亿 Windowing: Soft Tissue 🗸 360 60	- D × Phantom: BODY v
	D Image Ima
Sort Images < 1/48 > 1 Go to slice	

Figure 22. IndoseCT screen display for patient images and image processing results.

The <u>third part</u> is a group of buttons at the top of the image and at the bottom of the image (**Figure 23**). The upper group of buttons is for opening the image, displaying DICOM info, and specifying the image windowing to use. **Figure 24** shows a button to download the **IndoseCT** manual book. The bottom group of buttons is for deleting images, sorting images, and shifting images from one slice to another.

Figure 23. The third part displays a group of buttons at the top of the image and at the bottom of the image for various purposes, such as opening an image, displaying DICOM info, and so on.

G IndoseCT v20b		- 🗇 🗙
0 -	English Bahasa Indonesia ID Exam Date 16/09/2021 Name Institution Age - Scanner Manufacturer Sex Unspecified V Protocol	-Model
200 -	Methods: Calculation One slice	~
300 -	Parameter Output mAs Scanner Emotion Voltage (k/) 110 CTDL (mGy)	0.00
400 -	Tube Current (mA) 100.00 C1DL _{ad} (m6y) Rotation Time (s) 1.00 Pitch 1.00 Collimation (mm) 1	0.00
500 -	Scan Length (cm) 10.00 Calculate Get Info	Next
-100 0	100 200 300 400 500 600 701 0 Go to slice	

Figure 24. Button to download the IndoseCT manual book.

The <u>fourth part</u> for displaying patient information, such as ID, name, age, gender, examination time, hospital, type of scanner, and the protocol used (Figure 25). This information will be filled in automatically when the image is opened. If we calculate the dose without using the image and this information is to be stored in the database, then this information must be filled in manually.



Loading Images

Figure 25. Part to display patient information, such as ID, name, age, gender, examination time, hospital, type of scanner, and protocol used.

The <u>fifth part</u> is the graphical part. For every calculation, especially those involving profile data or data sets that have been stored in the database, the results can be displayed in a graph (Figure 26). These charts can be saved and the data can be exported to Microsoft Excel.



Figure 26. The Graphical display to view profile data of a parameter or data set that has been stored in the database. These charts can be opened and closed without interrupting calculations.

IV. PATIENT IMAGE

4.1. Open image

The advantage of **IndoseCT** is that $CTDI_{vol}$, D_{eff} , D_w , SSDE, and DLP can be taken and calculated from patient images. However, all these parameters can also be obtained without patient images. **IndoseCT** provides the option that all of these parameters can be obtained independently of the patient image. If we want to get the values of all these parameters from the patient image, then the patient image must be opened first.

IndoseCT has two options for opening images, namely by **file** and by **folder**. For educational purposes, **IndoseCT** has provided one sample, namely an anthropomorphic phantom from the lower end of the pelvis to the top of the head. The sample images can be opened by left-clicking the **Open Sample** button (**Figure 27**). In this sample there are 88 slices (images).



Figure 27. Open Sample button to open sample images.

To open an image by file, press the **Open File(s)** button at the upper left corner (**Figure 28**), and the **Open File** dialog will appear. Select the folder where the patient image files are stored. At that time, the file names will be empty. Then select **All Files**, under **Type Files**, which is located at the very bottom (see **Figure 29**). After that, the file name of the patient image will be displayed.

		- 0) ×
	/ing: Soft Tissue V 360 60	Phantom: B)DY ~
Open File(s)		ID Exam Date 16/09/2021 Name Institution	
100 -		CTDIvol Diameter SSDE Organ Analyze Methods: <td></td>	
		Calculation One slice	~
200 -		Parameter Output Manufacturer Siemens V MAS	0.00
300 -		Scanner Emotion Effective mAs Voltage (kV) 110 CTDI _w (mGy)).00).00
		Tube Current (mA) 100.00 CTDL _{vel} (mGy) C Rotation Time (s) 1.00 DLP (mGy-cm) C	.00
400 -		Pitch 1.00 Collimation (mm) 1 ~	
500 -		Scan Length (cm) 10.00 Calculate Get Info	
		N	ext
-100 0 100	200 300 400 500 600 70	0	
Sort Images < 0/0 > 0 Go to slice			

Open File(s)



Open Files				×				Pha	ntom: BOD
→ · ↑ 📙 « s	SIEMENS ANANG SURYAN > SE3	✓ Č							
ganize 🛪 New fol	lder	8	· • •	0	ID		Ex	am Date 16/09/2021	
	Name	Date modified	Type	^	Age -		Sc	anner Manufacture	r-Model
This PC	Acculmana dir.	09/02/2021 15:09	DIP Eile		Sex Unspecified		~ Pr	otocol	
3D Objects	IM0	13/02/2019 1:55	File		cronul at	0005			
Desktop		13/02/2019 1:55	File		CIDIVOI Diame	ter SSDE	Organ	Analyze	
Documents	1M2	13/02/2019 1:55	File		Methods:				
Downloads	IM3	13/02/2019 1:55	File		Calculation				
Music	IM4	13/02/2019 1:55	File		One slice				
Pictures	IM5	13/02/2019 1:55	File		one bite				
📲 Videos	IM6	13/02/2019 1:55	File						
Local Disk (C:)	IM7	13/02/2019 1:55	File		Parameter			Output	
DATA 1 (D:)	IM8	13/02/2019 1:55	File		Manufacturer	Siemens	\sim	mAs	0.0
DATA 2 (E:)	IM9	13/02/2019 1:55	File		Scappor	Emotion		Effective mAs	0.0
Maturali	< <	13/02/2019 1-55	File	>	Sedimen	Linddon		CTDL (mGv)	0.0
File	name	All Files			Voltage (kV)	110	~	CTDL (mGu)	0.0
1.10		All Files			Tube Current (mA) 1	00.00		0.0
		DICOM Files (*.d	cm)		Rotation Time (s)		1.00	DLP (mGy-cm)	0.0
					Pitch		1.00		
					Collimation (mm)	1	~		
					Scan Longth (cm)		10.00		
					Scan Lengur (cm)		10.00		
					Calculate	Get Info	C		
									Nex
-100	0 100 200	300 400	50	0 600 7	roc				

Figure 29. Display of the Open Files dialog.

Note that in general, patient images consist of more than one slice (file) (in some cases more than 1000 slices). If you choose only one file, then just select one file, click **Open** or directly double-click on the file. To select multiple patient image files, click the first file, then press the **Shift key** on the keyboard and hold it, then press the key to move the cursor down. Or the user can also press **Ctrl + A**. However, it must be remembered, that at the bottom of the file there is usually an **Acculmage.dir** file, this file should not be selected (**Figure 30**). When several files have been selected, the next step to open the image is to press the **Open** button. As noted the selected image is the patient's axial image file. Check that the axial image files

have the same file size, about 516 KB, and do not select files of a different size (These are not patient axial image files).



Figure 30. Display of selected patient image files (blue blocks).

Pressing the **Open** button will open all files. To open multiple images (of hundreds or thousands) usually takes a few seconds. On the computer screen, there is a notification of the progress of the file opening (**Figure 31**). Wait a moment for the axial images of the patient to be displayed on the computer screen.



Figure 31. Notification of patient axial images opening as a percentage of the total.

To open images by folder, select the **Open Folder** button which is to the right of the **Open Folder** button (**Figure 32**). When the button is pressed, the **Open Folder** dialog will appear. Select the folder where the patient image files are stored.



Figure 32. Open File(s) button to open an image by file.

Note that opening multiple images on a folder basis is more practical than on a file basis. However, it must be ensured that in the folder there are no other files except the patient's axial image file. If there are other files, then DICOM info from the axial image will not be able to be displayed later. This causes the next process (such as calculating patient diameter) to fail since calculating this parameter requires information from the DICOM info. An explanation of DICOM info is discussed in the next sub-chapter. An example of an opened patient image is shown in **Figure 33**.

🗅 🗅 📽 🅸 🕐 🐻 🖸 Wind	wing: Soft Tissue ~ 360 6	50					Phan	tom: BODY
			ID			Exam Da	te 28/01/2019	
			Name			Institutio	n	
			Age	39 M		Scanner	SIEMENS-Sen	sation 64
			Jex	IM		- FIOLOCOI	ABDOMEN	
			CTDI	vol Diamet	er SSDE	Organ Ar	alyze	
			Metho	ods:				
		5000	Calcu	ulation				
	and the same higher		One	slice				
			Para Mar Sca Volt Tub Rota	ameter nufacturer nner age (kV) ne Current (mA) ation Time (s)	Siemens Emotion 110 1	Outp mAs Effec CTD) 00.00 DLP	tive mAs	0.1 0.1 0.1 0.0
Contraction of the second	Same		Pitcl	h imation (mm)	1	1.00 ~		
			Sca	n Length (cm)		10.00		
				Calculate	Get Into			
								Nex

Loading Images

Figure 33. Example of patient image display that has been successfully accessed. The ID and name of patient and name of institution are erased for patient privacy.

4.2. Show image

Access to the next file (slice) is achieved by clicking the > and < buttons (below the image). The > button is for viewing the next slice, and the < button is for viewing the previous slice. In the middle is the number of slices visible on the screen and the total number of slices. In the example (**Figure 34**), the opened slice is number 30, and the total number of slices is 88.



Figure 34. The > button to move to the next slice, and the < button to move to the previous slice.

If the slice number is the same as the total number of slices, by clicking the > button will return the visible slice to slice number one. Conversely, if the slice number is equal to 1, then clicking the < button will display the last slice.

To quickly access a particular slice, you can fill in the number of the desired slice, and click **Go to slice**. This will greatly help speed up thedisplay, if we are working with hundreds or even thousands of images (slices), as is common with modern CT machines.

More practically, the user can use the arrow keys on the **keyboard**. The **up arrow key** is to view the next slice, the **down arrow key** is to view the previous slice, the **left arrow key** is to view the previous 5th slice, and the **right arrow key** is to view the next 5th slice.

In fact, often the images we get are not sequential, for example the first slice shows the head image, the second slice shows the chest image, the third slice shows the head image again, and so on. This will be very difficult. To **sort images**, you can do this by pressing the **Sort Images** button (at the bottom corner). Images can be deleted by pressing the **Close Images** button, which is to the left of the **Sort Images** button.

The CT number within images has a very wide range from about -1000 HU to over +3000 HU. With this wide range, if the image is displayed as it is, the contrast between objects will appear very low. To improve the contrast appearance of objects in image, the windowing technique is used (more information about windowing in books on CT). **IndoseCT** is equipped with **Windowing** selection, namely: **Soft Tissue**, **Bone**, **Lung**, **Liver**, **Brain**, **Fat**, **Spline**, **PF** (**Posterior Fossa**), **IAC** (**Internal Acoustic Canal**), **Vascular**, **Custom**, and **None**. The default is **Soft Tissue**. **Custom window** means that the user can define their own window-width (WW) and windowlevel (WL) values. **None** means that the image is displayed without using windowing. This means that the image is displayed from the minimum to the maximum value. The sample image when displayed using the **Bone** window looks like in **Figure 35** and using the **Lung** window looks like in **Figure 36**.

G IndoseCT v20b	- 0 ×
E D D T & C Windowing: Bone V 2000 400 Ph	antom: BODY ~
ID 999999-001 Exam Date 09/02/2017 Name XXXXXXX-001 Institution Kyushu.Uhh Age - Scance TOSHIBA-A Sex Unspecified Protocol -	exion
CIUVoi Diameter SSDE Organ Analyze Methods: Calculation One sice	~
200 - Parameter Output Manufacturer Siemens V MAS	0.00
Scanner Emotion CTL_(mGy) Tube Current (mA) 100.00 DUP (mGy) CTL_ue (mGy) Tube Current (mA) 100.00 DUP (mGy) CTL_ue (mGy) DUP (mGy) CTL_ue (mGy)	0.00 0.00 0.00 0.00
Kotation Irme (s) 1.00 Pitch 1.00 Collimation (rms) 1	
400 - Scan Length (cm) 10.00	
S00 Calculate Get Into -100 0 100 200 300 400 500 600	Next

Figure 35. Sample image of the chest area displayed with the Bone window.

					- 0 ×
				1	-nancom. BODT V
	ID 999999-001		Exa	am Date 09/02/20	17
0 -	Name XXXXXX-001		Ins	titution Kyushu.U	Alexion
	Sex Unspecified		V Pro	itocol	Alexion
	CTDIvol Diamet	er SSDE	Organ	Analyze	
100	Methods:				
	Calculation				~
	One slice				~
200 -	Parameter			Output	
	i di difietei			mAs	0.00
	Manufacturer	Siemens	~	Effortivo mác	0.00
	Scanner	Emotion	~	CTDL (mos)	0.00
200	Voltage (kV)	110	~	CTDI _w (mGy)	0.00
	Tube Current (mA)	100.00	CIDI _{vel} (mGy)	0.00
	Rotation Time (s)		1.00	DLP (mGy-cm)	0.00
	Pitch		1.00		
	Collimation (mm)	1	~		
400 -	Scan Length (cm)		10.00		
	Sour Lenger (em)		10.00		
	Calculate	Get Inf	fo		
500					
					Next
-100 0 100 200 300 400 500 600	1.00				
Sort Images < 30/88 > 30 Go to slice					

Figure 36. Sample image of the chest area displayed with the Lung window.

4.3. Zoom-in, zoom-out and move image

The opened patient image can be enlarged (zoom-in) or minimized (zoom-out) or shifted as needed. To enlarge or minimize an image, place the cursor over the image, then use three fingers to enlarge the image. An enlarged anthropomorphic phantom image is shown in **Figure 37** and a minimized image is shown in **Figure 38**.

10 10 99999001 Exam base 90/02/017 10 Name 2000000000000000000000000000000000000	indoseCT v20b	– Phantom: [日 > BODY
		D 999999-001 Exam Date 09/02/2017 Name XXXXXX-C001 Institution Kyushu.Univ Age	0.00 0.00 0.00 0.00



G Ind	oseCT	v20b		•															-	. 0	×
لغا			4 XX	: (?)		≌	Windowing	: Soft Tissu	e ~ 360	60									Phan	tom: BODY	~
c	_													ID Name Age Sex	999999-001 XXXX-001 - Unspecified		v DE Orr	Exam Date Institution Scanner Protocol	09/02/2017 Kyushu.Univ TOSHIBA-Alex	ion	
100														Metho Calcu One s	Ids: Idation	er ssi	JE OR	jan Ana	yze		~
200 300								6						Para Man Scar Volti	ameter ufacturer nner age (kV)	Siemens Emotion 110	s ~ ~	Output mAs Effectiv CTDI _w	re mAs	0.00	
400														Tub Rota Pitch	e Current (mA) ation Time (s) n		100.00 1.00 1.00	DLP (r	nGy-cm)	0.00]
500														Colli Scar	mation (mm) n Length (cm) Calculate	1 Get					
600		00	-10	0	0		100	200	300	400	500	600	700							Next	
	Sort	Images	<	30/88	> 30) Go	to slice														

Figure 38. Image of a minimized sample of the chest area.

To move an image, place the cursor over the image, hold down the left-click and shift the image to the desired position. An example of a shifted image is shown in **Figure 39**.

G IndoseCT v20b		- o >
┣┗┗뿍‡\$?. ВС	Windowing: Soft Tissue V 360 60	Phantom: BODY
100 -		ID 999999-001 Exam Date 09/02/2017 Name XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
		CTDIvol Diameter SSDE Organ Analyze
		Methods:
200		Calculation ~
		One slice ~
300 -		Parameter Output Manufacturer Siemens mAs 0.00 Scanner Emotion Effective mAs 0.00 Voltage (kV) 110 CTDI _w (mGy) 0.00
400		Tube Current (mA) 100.00 CTDI _{vel} (mGy) 0.00
400 -		Rotation Time (s) 1.00
		Pitch 1.00
		Collimation (mm) 1 ~
500		Scan Length (cm) 10.00
		Calculate Get Info
600 -		Net
-300 -200	-100 0 100 200 300 400	Next
Sort Images < 30/88 > 30 G	o to slice	-

Figure 39. Sample image of the displaced chest area.

4.4. Image pixel value

When taking D_{eff} or D_w measurements, sometimes we want to know the pixel value at a certain point in the patient's image. In IndoseCT, to find out the pixel value and its position (at the x- and y-coordinates), simply place the cursor on the pixel in question. The pixel value and its position will be displayed dynamically at the top of the image, as shown in Figure 40.



Figure 40. When the cursor position is placed at a certain position in the image, the pixel value will be obtained and displayed automatically at the top of the image.

4.5. DICOM info

Unlike non-medical images, medical images are saved in the DICOM format. In the DICOM format, the image also contains detailed information about the image, for example image size, number of pixels, date of data acquisition, patient name, patient age, and so on. This information is often referred to as DICOM info or DICOM header. The DICOM details may vary from one scanner to another. To display DICOM info, press the **DICOM Info** button. An example of DICOM info can be seen in **Figure 41**. To find out the meaning of DICOM info, refer to books or journals that discuss DICOM info. The user may need to retrieve some of this information for certain purposes, for example if the user wants to know the height of the patient table to calculate mis-centering or the user wants to know the field of view (FOV) of the image obtained.

Note that the DICOM info that appears is for the currently active image on the screen. If we want to see DICOM info for the next slice, then the image must be moved to the next slice first. In general, the DICOM info for each slice is similar except for some parameters such as slice location, tube current value (for TCM applications), and image acquisition time.

		ID 999999-001 Exam Date	09/02/2017
		Name XXXXXX-001 Institution	Kyushu.Univ
	G DICOM Tree View	? ×	TOSHIBA-Alexion
	1	^	
	(0008, 0008) Image Type	CS: ['ORIGINAL', 'PRIMARY', 'AXIAL']	rze
	(0008, 0016) SOP Class UID	UI: CT Image Storage	
	(0008, 0018) SOP Instance UID	UI: 1.2.392.200036.9116.2.5.1.37.2427178946.1486599611.328570	
	(0008, 0020) Study Date	DA: '20170209'	
	(0008, 0021) Series Date	DA: '20170209'	
	(0008, 0022) Acquisition Date	DA: '20170209'	
	(0008, 0023) Content Date	DA: '20170209'	
	(0008, 0030) Study Time	TM: '090946.000'	
	(0008, 0031) Series Time	TM: '091346.669'	
	(0008, 0032) Acquisition Time	TM: '091401.850'	
	(0008, 0033) Content Time	TM: '091415.078'	
	(0008, 0050) Accession Number	SH: '589'	e mas
	(0008, 0060) Modality	CS: 'CT'	mGy)
	(0008, 0070) Manufacturer	LO: 'TOSHIBA'	(200)
	(0008, 0080) Institution Name	LO: 'Kyushu.Univ'	(moy)
	(0008, 0090) Referring Physician's Name	PN: ''	Gy-cm)
	(0008, 1010) Station Name	SH: 'ID STATION'	
	(0008, 103e) Series Description	LO: ',, Axial, 7.0, '	
	(0008, 1040) Institutional Department Name	LO: 'ID DEPARTMENT'	
	(0008, 1090) Manufacturer's Model Name	LO: 'Alexion'	
	(0010, 0010) Patient's Name	PN: 'XXXXXX-001'	
	(0010, 0020) Patient ID	LO: '999999-001'	

Figure 41. Button to display DICOM info and an example of DICOM info.

4.6 Patient data

Before calculating patient's dose, patient data be filled in first. This data is useful for analysis at the final stage. The data that needs to be filled in includes (see **Figure 42**):

- Patient ID (ID)
- Patient name (Name)
- Patient age (Age)
- Patient's gender (Sex)
- Date of examination (Exam Date)
- Institution where examination is performed (name of hospital) (Institution)
- Scanner type, both manufacturer and scanner model (Scanner).

• Examination protocol (**Protocol**)





Note: When the patient image is opened, the data input parameters above are directly extracted from the DICOM info of the image. An example of patient information that has been filled in when opening a sample image can be seen in **Figure 43**.



Figure 43. Filled form from DICOM info of the sample image. **Age**, **Sex** and **Protocol** content fields are not filled because the information was not available in the DICOM info of the sample image.
V. DETERMINATION OF THE CTDIvol

The magnitude of the CT dose is generally expressed using the metric $CTDI_{vol}$. $CTDI_{vol}$ is a quantity that determines the magnitude of the output dose of the CT machine. The $CTDI_{vol}$ value is strongly influenced by exposure factors such as tube voltage (kVp), tube current (mA), rotation time (s), pitch, and collimation width. The $CTDI_{vol}$ is specific for each scanner. The patient dose is calculated based on the magnitude of the $CTDI_{vol}$ combined with the specific characteristics of the patient.

The second quantity that describes the total energy transferred to the patient is the doselength product (DLP). This DLP is calculated based on the product of CTDI_{vol} and scan length (L). Calculation of the patient's effective dose is usually calculated from this DLP and conversion factor (f).

On modern CT machines, these CTDI_{vol} and DLP values are usually displayed on the CT console screen. For the newest product CT machines, these two values are also stored in DICOM info. In addition, the CTDI_{vol} and DLP values are also quantities that must be checked periodically in the quality control (QC) program.

To calculate $CTDI_{vol}$, first select the $CTDI_{vol}$ tab in the row of tabs. If the $CTDI_{vol}$ tab is clicked, the color will change from gray to white. By default, when IndoseCT is activated, this $CTDI_{vol}$ calculation appears on the first screen.

In IndoseCT 20.b, there are three options for determining it: calculation (Calculation), manual (Input Manually), and taken from DICOM info (Get From DICOM) (Figure 44). The default option is Calculation.



Figure 44. There are three options to obtain the CTDI_{vol} value in IndoseCT, namely Calculation, Input Manually, and Get from DICOM.

5.1 Input manually

When the user already has the $CTDI_{vol}$ and DLP values, either from the CT console screen or from measurements (e.g. when conducting a quality control (QC) program or compliance test) or from other software, the user enters these values into this **IndoseCT**. The steps are as follows (see **Figure 45**):

- Select Enter Manually option.
- Fill in the **CTDI**_{vol} (mGy) box with the CTDI_{vol} value.
- Fill in the **DLP (mGy-cm)** box with the DLP value.

G indoseCTv20b	-	٥	\times
🖿 🖿 🥵 🏟 🕐 🖬 🔀 Windowing: Soft Tissue ~ 360 60	Phantom:	BODY	~
0 - D D 99999-001 Exam Date (0/02/ Name XXXXX-001 Institution Kyushu Age - Scanner TOSHI Sex Unspecified V Protocol	/2017 J.Univ IBA-Alexion		
100 - CTDIvol Diameter SSDE Organ Analyze Methods: Input Manually		~	
200 - CTDL _{ref} (mGy)		0	1
300 Submit			
400 -			
		Next	
Start Images < 30/88			



5.2 Calculation

a. One slice

The default $CTDI_{vol}$ and DLP calculations on **IndoseCT** is **One slice** option (Figure 46). In addition to the **One slice** option, another option is **Z-axis**.

If the user does not yet have the CTDI_{vol} and DLP values, then the user can calculate these two quantities by entering the type of CT used for the exposure factors. In **IndoseCT**, the only available options are for CT manufacturers **Siemens**, **Toshiba**, **General Electric** (GE) and **Philips** (See **Figure 46**). For other manufacturers, values are still in process.

C₀ IndozeCT v20b	- 0" × Phantom: BODY v
0 -	ID Exam Date 16/09/2021 Name Institution Institution Age Unspecified V Protocol CTDIvol Diameter SSDE Organ Analyze Methods: Calculation V
200 - 300 -	One sice Output Parameter Output Manufacturer Siemens Scanner Toshiba Voltage (kV) GE Tube Current (mA) 100-00 Ubit Current (mA) 100-00
400 -	Rotation Time (s) 1.00 Pitch 1.00 Collimation (mm) 1 Scan Length (cm) 10.00
-100 0 100 200 300 400	Calculate Get Info Next

Figure 46. There are four manufacturer options in IndoseCT.

The steps to calculate CTDI_{vol} are as follows:

- Select the **CTDI**_{vol} tab (Actually, when **IndoseCT** is opened, the first tab is **CTDI**_{vol}).
- Select the **Calculation** option (Actually, when IndoseCT is opened, the first option is to specify this CTDI_{vol}).
- Select the Manufacturers name.
 There are four options including Siemens, Toshiba, General Electric (GE) and Philips.
 The default is Siemens.
- Select the type of **Scanner**.
 - There are 11 types of scanners for Siemens: Emotion, Emotion Duo, Emotion 6, Sensation 4, Sensation 10, Sensation 16, Sensation 16 Straton, Sensation 64, Sensation Open, Somatom Definition AS, and Somatom Perspective. A selection of Siemens scanner types is shown in Figure 47.
 - There are 8 scanners for Toshiba: Asteion, Asteion Dual, Asteion Multi, Aquilion Multi/4, Aquilon 16, Aquilion 64, Alexion, and Alexion Access.
 - There are 7 scanners for GE: HiSpeed, HiSpeed ZX/i, LightSpeed, LightSpeed Plus, LightSpeed Ultra, LightSpeed 16 and LightSpeed VCT.
 - There are 9 types of scanners for Philips: AV, LX, M/EG, Mx8000, SR700, CT Secura, Brilliance 16, Brilliance 64, and Ingenuity.

The default scanner is Emotion Duo.

- Select Voltage (kV). The voltage value depends on the type of scanner used.
- Fill in the **Tube Current (mA)** value. The default tube current value is **100 mA**.
- Fill in the **Rotation Time(s)**. The default rotation time is **1 s**.

	– 0 ×
Windowing: Soft Tissue V 360 60	Phantom: BODY ~
0 - Dameter SSDE Organ	am Date 16/09/2021 ttution manufacturer-Model Analyze
200 -	Output mAs 0.00 Effective mAs 0.00
300 - Voltage (k/) Tube Current (mA) Rotation Time (s) Sensation 16	CTDI _{vel} (mGy) 0.00 CTDI _{vel} (mGy) 0.00 DLP (mGy-cm) 0.00
Sensation 10 Sensation 10 Sensation 40 Sensation 40 Sensation 40 Sensation 64 Sensation 64	
-100 0 100 200 300 400 500 600 700	Next

Figure 47. A list of scanner options available for Siemens on IndoseCT.

- Fill in the **Pitch**. The default pitch is **1**.
- Select **Collimation Width (mm)**. The choice of collimation width is determined by the manufacturer's name and the type of scanner.
- Fill in the Scan length (cm). The default scan length is 10 cm.
- When everything is filled in, then press the **Calculate** button.

The CTDI_{vol} and DLP values will then be calculated. In addition, the mAs value, effective mAs (mAs divided by pitch), and CTDIw values will also be obtained. The CTDI_{vol} value is CTDI_w divided by pitch, and DLP. An example of CTDI_{vol} and DLP calculation results for Siemens Emotion is shown in **Figure 48**.

For practical purposes, input parameters for calculating CTDI_{vol}, such as tube current, tube voltage and so on, can also be taken from the DICOM info. However, this only applies if the CT manufacturer and scanner type are already in the existing scanner database. To retrieve DICOM info, press the **Get Info** button. For example, the input parameters taken from the DICOM info for the Siemens Sensation 64 are shown in **Figure 49**.

G IndoseCT v20b	– ō ×
🖿 🖿 🦉 🤻 🛞 💽 🔀 Windowing: Soft Tissue 🖌 360 60	Phantom: BODY ~
0 - ID Exam D Name Institut Age - Scame Sex Unspecified ~ Protoco	Date 16/09/2021 Jtion
100 - CTDIvol Diameter SSDE Organ A Methods:	Analyze
200 -	~
Parameter Out Manufacturer Siemens V	itput iAs 100.00
Scanner Emotion Voltage (kV) 110 V	rective mAs 100.00 FDI _w (mGy) 0.00 TDI_{vol} (mGy) 0.00
400 - Pitch 1.00	LP (mGy-cm) 0.00
Collimation (mm) 1 v Scan Length (cm) 10.00	
500 Calculate Get Info	
-100 0 100 200 300 400 500 600 700	Next
Sort Images < 0/0 > 0 Go to slice	

Figure 48. CTDI_{vol} and DLP calculation results for Siemens Emotion.

h IndoseCT v20b Im Im I	- 🗗 🗙 Phantom: BODY 🗸
	ID Exam Date 28/01/2019 Name Institution Institution Sex V Protocol Sex V Protocol CTDIvol Diameter SSDE Organ Analyze Methods: Image: Calculation Calculation Voltage Voltage Senaner Senaner Senaner Senaner Senaner Voltage (N) Nucho Current (mA) 178.00 Pitch 1.400 Collimation (mm) 28.8 Scan Length (cm) 11.00 Calculate Get Info
100 0 100 200 300 100 000 Images < 1/48 > 1 Go to slice 000 000 000	

Figure 49. Input parameter values are retrieved from DICOM Info with the **Get Info** button to calculate CTDI_{vol}.

For CT applications using the tube current modulation (TCM) technique, the **Tube Current** (mA) is not known at the beginning of scanning. The value of tube current (mA) varies greatly along the z-axis and is determined during scanning based on the patient's scout image (scanogram or localizer). If a patient image is available, the **Tube Current (mA)** value can be extracted from the DICOM info for the patient image. This **IndoseCT** software has also been equipped with a feature to extract these varied mA values, namely by pressing the **Get Info** button.

b. Z-axis option

In order to display the tube current rating (mA) info and calculate the CTDI_{vol} value along the longitudinal-axis, the **Z-axis mode** option must be selected **(Figure 50**).

When **Z-axis mode** is selected there are three options, namely **Slice Step**, **Slice Number**, and **Regional (Figure 50)**. If **Slice Step** is selected and the field is filled with 1, it means that the CTDI_{vol} value is calculated from all slices. If the field is filled with 2, it means that the CTDI_{vol} value is calculated from the 1st slice, then the 3rd slice, then the 5th slice, and so on.

If **Slice Number** is selected and the field is filled with 1, it means that the $CTDI_{vol}$ value is calculated from the (single) middle slice. If the field is filled with 3, it means that the $CTDI_{vol}$ value is calculated from one slice in the middle, one slice at position 1/3 from the beginning, and one slice at position 1/3 from the end.

Meanwhile, if **Regional** is selected, then we can specify a slice to start (e.g. the 15th slice) and a slice to end (e.g. the 25th slice). **Regional** selection is very useful for estimating the dose of organs located at certain positions. It has been reported from several studies that organ dose estimation using CTDI_{vol} calculated from the regional mean is more accurate than that calculated from the average of all existing slices (Anam et al. J Biomed Phys Eng. 2021; Khatonabadi et al. Med Phys. 2013; 40(5): 051905; Bostani et al. Med Phys. 2015; 42(2): 958-968).



Figure 50. Z-axis mode for calculating tube current values (mA) and $CTDI_{vol}$ along the longitudinal axis for the TCM technique. Tube current rating (mA) is taken from the DICOM info.

The tube current profile graph and CTDI_{vol} can be displayed by selecting **Show CTDI_{vol} Graph** and **Show mA Graph** (**Figure 50**). The two graphs can be displayed simultaneously or separately. After all the options are determined, the **Calculate** button is pressed, and the tube current (mA) (**Figure 51**) and CTDI_{vol} values will automatically be displayed along the

longitudinal-axis (Figure 52). The average value is filled in the Tube Current (mA) and CTDI_{vol} (mGy) fields.



Figure 51. Graph of tube current profile (mA) along the longitudinal axis for the TCM technique. This tube current (mA) is extracted from the DICOM info of patient image.



Figure 52. Graph of the $CTDI_{vol}$ (mGy) profile along the longitudinal axis for the TCM technique.

The average tube current (mA) value and its standard deviation can be displayed on the graph by pressing the **Option** button. Then proceed by selecting the **Mean** and **Standard Deviation**

of the y-data value. The $CTDI_{vol}$ (mGy) profile graph with the mean and standard deviation is shown in **Figure 53**.



Figure 53. Graph of the CTDIvol (mGy) profile along the longitudinal axis for the TCM technique with the mean and standard deviation lines.

The x-axis and y-axis of the profile graph can be set by the user to the minimum and maximum values. Right-click on the graph, then highlight the x-axis or y-axis, then the minimum and maximum values can be determined manually to make the graph look optimal. The tube current graph (mA) for which the minimum value of the y-axis has been set to **100 mA** is shown in **Figure 54**.

The resulting graph can be saved (Save Plot) in several formats. Write file name, and select Save as type to choose PNG, TIFF, JPEG, Bitmap, Scalable, Comma-Separated Value, or Microsoft Excel Workbook (Figure 55).



Figure 54. Graph of tube current (mA) where the minimum value of the y-axis has been changed to **100 mA**.



Figure 55. Options for saving tube current profiles and $CTDI_{vol}$ along the longitudinal, by selecting Save as type.

An example of a Microsoft Excel data/file from a $CTDI_{vol}$ profile along the longitudinal-axis is shown in **Figure 56**. We can then manipulate the data in Microsoft Excel as needed.

- 🗇 🗧 = IndoseCT - 🗇 👌	×
File Home Share View	• 🕐
← → * ↑ → This PC > DATA 2(E) > Pasca ITB > IndoseCT ♥ Ŏ Ø Search IndoseCT	
 	
12 items	

Figure 56. Example Microsoft Excel file of a CTDI_{vol} profile along the longitudinal axis.

5.3 From DICOM info

For the latest CT devices, usually $CTDI_{vol}$ is already stored in DICOM info. In this case, the $CTDI_{vol}$ value can be directly taken from the DICOM info, and the DLP value is calculated as the product of the $CTDI_{vol}$ value and the scan length (L). The scan length itself is calculated as the distance between the slice location for the last slice minus the slice location for the first slice, then add the beam width.

a. One slice

To extract $CTDI_{vol}$ from the DICOM info select **Get from DICOM** from the options. Choose whether only one slice (**One slice**) or **Z-axis**.

For **One slice**, press the **Calculate** button (**Figure 57**), and the CTDI_{vol} value will be obtained for that slice (**Figure 58**).

If the $CTDI_{vol}$ value is already available in the DICOM info, **Get from DICOM** option is very straightforward. However, in reality, not all DICOM info stores $CTDI_{vol}$ values. If the DICOM info does not store the $CTDI_{vol}$ value, then **IndoseCT** will assign a value of 0 to $CTDI_{vol}$ (**Figure 58**). Thus, we must use one of the other two options to obtain the $CTDI_{vol}$ value.

Windowing 🗈 🗅 📽 🏟 🕐 🗟 🖸	: Soft Tissue ~ 360 60			Phantom: BODY
		ID		Exam Date 28/01/2019
0		Age 3	39	Scanner STEMENS-Sensation 64
		Sex M	м ~	Protocol ABDOMEN
		СТДІУО	Diameter SSDE Org	an Analyze
100		Get from	s: m DICOM	v
	and the second se	One slic	ce	~
200 - 300 -		Param Scan L CTDI _u DLP (m Calcu	heter Length (cm)	10] 0] 10 10 10
500				Next

Sort Images

Figure 57. The option selected for calculation on the visible slice is **One slice**. After that proceed by pressing the **Calculate** button.

C. IndosecT v20b	0		- D X Phantom: BODY V
0-		ID Name Age 39 Sex M CTDIvol Diameter SSDE O Methods: Get from DICOM One slice	Exam Date 28/01/2019 Institution Scanner SIEMENS-Sensation 64 ABDOMEN rgan Analyze
300		Parameter Scan Length (cm) CTDI _{vd} (mGy) DLP (mGy-cm) Calculate Adjust CTDIvol	11.00 4.29 47.16 with mA
400 500 -100 501 Images < 1/48 > 1 Go to slice Locium Images	400 500 600		Next

Figure 58. For One slice mode, the $CTDI_{vol}$ value in the slice is obtained.

Figure 59. If DICOM info does not store the $CTDI_{vol}$ value, then **IndoseCT** will assign a value of 0 to $CTDI_{vol}$ and DLP.

In the TCM technique, the tube current (mA) varies for each slice, causing the $CTDI_{vol}$ value to vary also. In reality, not all $CTDI_{vol}$ values when using the TCM technique will vary according to the tube current value. For this situation, the $CTDI_{vol}$ value can be adjusted by selecting **Adjust CTDI**_{vol} with mA (Figure 60).

Image: Serie Piet I is a Piet I is a Piet I issue is and is an interview is a construction of the serie is construction of the serie is a constr	G IndoseCT v20b			- 0 ×
0 100 200 <	┣┗ ┗ ╬₡Ø. ∎	Windowing: Soft Tissue v 360 60		Phantom: BODY ~
	0 - 100 - 200 - 300 - 500 -		ID Exam Date 28/01 Name Institution Age 39 Scanner SIEM Sex M ✓ Protocol ABDC CTDIvol Diameter SSDE Organ Analyze Methods: Get from DICOM One slice Parameter Scan Length (cm) CTDL _w (mGy) DIP (mGy-cm) Calculate ✓ Adjust CTD2vol with mA	/2019 ENS-Sensation 64 MEN

Figure 60. Not every CT machine that uses the TCM technique produces a CTDI_{vol} value that follows the tube current value. For this situation, the CTDI_{vol} value can be adjusted by selecting **Adjust CTDI_{vol} with mA.**

b. Z-axis option

To find out whether the CTDI_{vol} that applies to TCM has been adapted or not, select the **Z-axis** option. There are then three further options, namely **Slice Step**, **Slice Number**, and **Regional** (**Figure 61**). **Figure 61** shows that the CTDI_{vol} value on the **Siemens Somatom 64** already follows the current value (mA).



Figure 61. The CTDI_{vol} value on the **Siemens Somatom 64** has been adjusted to the current value (mA).

Figure 62 shows that the $CTDI_{vol}$ value on the Toshiba Aquilio 128 does not follow the tube current value (mA), even though this CT machine uses the TCM technique. In **Figure 62**, the $CTDI_{vol}$ must be adjusted to the tube current value, namely by selecting **Adjust CTDI**_{vol} with **mA**. When this is done, the $CTDI_{vol}$ value automatically fluctuates according to the tube current value (mA) (**Figure 63**).

_				1	
	Graph of Slice - CTDIvol		- 🗆 ×	ID	Exam Date 11/11/2014
0		Sice - CTDIvol		Age 67	Scanner TOSHIBA-Aquilion
	°			Sex F	V Protocol CHEST
	5			CTDIvol Diameter SSDE	Organ Analyze
100	G Graph of Slice - Tube Current		– 🗆 X	Methods:	
		Chan Tuba Cumant		Get from DICOM	~
	َدُورِ اللهِ ال	Sice - Tube Current		Z-axis	~
00 -				Parameter	Z-axis Options
				Scan Length (cm)	37.60 Slice Step
	100			Avg. CTDI _{vol} (mGy)	5.70 O Slice Number
	(M)			DLP (mGy-cm)	214.32 O Regional
00	ant (n			Calculate Adjust CTDI	/ol with mA
	50 60				Show CI DIVol Graph
	Tube				
00	40				
	20				
	0				
00	0	100 200	300		Next
-100	0	Slice			

Figure 62. The CTDI_{vol} value on the Toshiba Aquilio 128 has not been pre-adjusted to the tube current (mA) value.



Figure 63. By selecting **Adjust CTDI**_{vol} **with mA**, the CTDI_{vol} value in the use of TCM which initially did not follow the tube current (mA) now fluctuates according to the tube current value (mA).

If the $CTDI_{vol}$ and DLP values have been filled in (using one of the three methods described, namely: **Calculation**, **Get from DICOM**, and **Input Manually**), then the user can proceed to the second stage, which is calculation of the effective diameter (D_{eff}) or water-equivalent diameter (D_{w}).

Note: The $CTDI_{vol}$ values obtained are based on two phantom sizes, i.e. either a phantom body (with a diameter of 32 cm) or a head (with a diameter of 16 cm). When determining the $CTDI_{vol}$ and DLP values, we must also carefully select the phantom type (phantom size) by selecting the **Phantom** option in the upper right corner. For $CTDI_{vol}$ values obtained automatically on the image, the phantom type has also been filled in automatically. However, the user needs to do double-check to ensure this selection. Incorrect selection of this type of phantom can result in a very large error in the estimated dose value.

VI. CALCULATION OF EFFECTIVE DIAMETER (D_{eff})

In principle, there are two types of diameter used to calculate the radiation dose received by the patient, namely the effective diameter (D_{eff}) and the water-equivalent diameter (D_w). The earlier D_{eff} is used for patient dose estimation, while D_w was introduced later and is more accurate. D_{eff} only takes into account the patient's geometry, while D_w takes into account the geometry and attenuation of the organ (organ composition) simultaneously.

IndoseCT can be used to calculate both D_{eff} and D_w . Select the **Diameter** tab, and choose either D_{eff} or D_w , as shown in **Figure 64**. (In this chapter, we will focus on discussing how to determine D_{eff} , then in the next chapter we will discuss on how to determine D_w).



Figure 64. To calculate the diameter, select the **Diameter** tab. There are two types of diameter that can be used, namely the effective diameter (D_{eff}) and the water equivalent diameter (D_w) .

To calculate the D_{eff}, there are two main options, namely the value is entered manually (**Input Manually**) or calculated from the image (**Get from Image**) as shown in **Figure 65**.

If the D_{eff} value is entered manually, the user should already have some patient data, such as AP diameter, LAT, AP+LAT, or patient age. Otherwise, if the D_{eff} value is taken from the image, the patient image in DICOM format must have been opened (how to open it has been discussed in the previous chapter).

	- 0 ×
Windowing: Soft Tissue V 360 60	Phantom: HEAD ~
0 - D Baneter SSDE Organ Ana Based on: Option	e 08/08/2021
200 - 300 - 300 -	hold (HU) 300 \$ txel area 500 \$ nod Lenter Aax
400 -	
500	Next
Soft Imagesi < 0/0 > 0 Go to slice	
READY	

Figure 65. There are two main options for calculating D_{eff} , namely the value is entered manually (**Input Manually**) or calculated from the image (**Get from Image**).

6.1 Manual calculation

D_{eff} value can be calculated manually by selecting the menu option of **Input Manually**. There are five options for **Input Manually** (Figure 66).



Figure 65. There are five options for calculating D_{eff} manually, using the diameter of **D_{eff}**, **AP**, **LAT**, **LAT+AP**, and **AGE**.

a. Input D_{eff} manually

If the user already has a D_{eff} value, then the D_{eff} value can be directly inputted to **IndoseCT**. In this method, no calculation is actually carried out. The steps are as follows (see **Figure 67**)

- Select Input Manually option.
- Select **D**_{eff} option.
- Fill in the value for the patient
- Press the **Calculate** button.

The D_{eff} value will be displayed in the D_{eff} (cm) box and visualized in a graph.



Figure 67. The display for inputting the D_{eff} value and the result. D_{eff} graph is shown with red dot.

b. D_{eff} calculation from AP diameter

To calculate the effective diameter (D_{eff}) from the AP diameter (see **Figure 68**):

- Select Input Manually option.
- Select **AP** option.
- Fill in the patient's AP diameter. (In this case, set the diameter value is set at 10 cm).
- Click Calculate button

The D_{eff} value will be calculated and the result will be displayed in the D_{eff} (cm) box, which is 11.60 cm. In this case, it is also displayed visually with graphics. The x-axis shows AP with a size of 10 cm, and the y-axis shows a D_{eff} with a size of 11.60 cm.



Figure 68. The display of the D_{eff} calculation from the AP diameter. The D_{eff} value in graph is shown by a red dot.

c. D_{eff} calculation from LAT diameter

Calculate the D_{eff} of the LAT diameter as follows (see Figure 69):

- Select Input Manually option.
- Select LAT option.
- Fill in the patient's LAT diameter. In this example, it is filled with 10 cm.
- Press the **Calculate** button.



Figure 69. The display for the D_{eff} calculation from the LAT diameter. The red dot in a graph shows the D_{eff} value.

After that, the D_{eff} value is calculated, and the results are displayed in the D_{eff} (cm) box, which is 10.20 cm. In this case, it is also displayed visually with graphics. The x-axis shows LAT with a magnitude of 10 cm, and the y-axis shows a D_{eff} with a size of 10.20 cm.

d. D_{eff} calculation from AP+LAT diameter

Calculate the D_{eff} from the AP+LAT diameter as follows (see Figure 70):

- Select Input Manually option.
- Select **AP+LAT** option.
- Fill in the patient's AP diameter. In this example, it is filled with 10 cm.
- Fill in the patient's LAT diameter. In this example, it is filled with 15 cm.
- Press the **Calculate** button.

After that, the D_{eff} value is calculated, and the results are displayed in the D_{eff} (cm) box, which is 12.20 cm. In this case it is also displayed visually with graphics. The x-axis shows AP+LAT with a magnitude of 25 cm, and the y-axis shows a D_{eff} with a size of 12.20 cm.



Figure 70. Display for D_{eff} calculation from AP+LAT diameter. The red dot in a graph indicates the D_{eff} value.

e. D_{eff} calculation from age

If the diameter of the patient is unknown, then the D_{eff} value can be estimated from the patient's age. However, it should be noted that the D_{eff} value based on age has a large variability. D_{eff} calculations with patient age should be avoided if possible, because its accuracy is low. Obtaining D_{eff} values from other methods is preferable.

Steps to estimate D_{eff} from the patient's age are as follows (see Figure 71):

- Select Input Manually option.
- Select AGE (max 18) option. In this example, it will be 5 years and 2 months.

• Fill in the patient's age in the **Year** box. Note: Maximum age is 18 years, and minimum is 0 years.



• Press the **Calculate** button.

Figure 71. Display for the D_{eff} calculation of the patient's age. The red dot in a graph indicates the D_{eff} value.

After that, the D_{eff} value is calculated, and the results are displayed in the D_{eff} (cm) box, which is 18.57 cm. In this case, it is also displayed visually with graphics. The x-axis shows the patient's age 5.17 years (5 years 2 months), and the y-axis shows D_{eff} with a size of 18.57 cm.

Note: When calculating the D_{eff} value using age, the **Phantom** must be **the Body Phantom** (not the **Head Phantom**)

6.2 D_{eff} calculation from image

In this case, the D_{eff} value is carried out directly from the patient's axial image. The user must select **Get from Image** from the available options, after ensuring that the patient's CT image has been opened. There are three choices of methods, namely: **Manual**, **Auto**, and **Auto** (**Z**-**axis**) (Figure 72). We will start with the **Manual** method.

0 Image Ima) 🗅 🗅 📲 🏟 🕐 💽 🚺 Windowing: Soft Tissue 🗸 360 60		Phantom: BODY
0 - - Institution 100 - - - - 200 - - - - 200 - - - - 200 - - - - 200 - - - - 200 - - - - 200 - - - - - 200 - - - - - - 200 - <		ID	Exam Date 28/01/2019
100 Image: Section of Ambroad State (State (Sta		Name	Institution
10 10 10 20 10 20 10 20 10 20 10 10 20 10 <td< td=""><td></td><td>Age 39 Sex M</td><td>Protocol ABDOMEN</td></td<>		Age 39 Sex M	Protocol ABDOMEN
100 200 300 400 -100 0 100 200 300 400 500 400 500 500 500 500 5		errer i Diemeter opprige	
100 200 300 400 -100		CTDIvol Diameter SSDE Or	jan Analyze
200 300 400 	0	Based on:	Options
200 - 300 - 500 - -100 0 100 200 300 400 500 600		Effective Diameter (Deff) \sim	AP 0 cm
200 300 400 		Source:	LAT 0 cm
300 -<		Get from Image ~	Clear
300 Image: Constraint of the second		Manual	1
300 -<		Auto	
300 300 300 300 300 300 400 500 600		Auto (Z-axis) Manual	
000 0 100 200 300 400 500 600		Show Graph	4
000			
000			
000 0 100 200 300 400 500 600 Previous Next			
000 -100 0 100 200 300 400 500 600			
-100 0 100 200 300 400 500 600			
-100 0 100 200 300 400 500 600			
-100 0 100 200 300 400 500 600			
-100 0 100 200 300 400 500 600		Previous	Next
	-100 0 100 200 300 400 500 60	00	

Figure 72. D_{eff} calculation from the image. There are three choices of methods, namely: **Manual**, **Auto**, and **Auto (Z-axis)**. When the user uses the **Get From Image** option menu, the patient's CT image in DICOM format must have already been opened.

a. Manual calculation

To perform a manual D_{eff} calculation from an image, steps are as follows (see **Figure 73**):

- First select **Manual** from the available options.
- Next, two buttons and two fields/in box will appear.
- Press the LAT button, a line will appear on the patient image in the lateral direction. We have to shift the line and make sure the left edge is on the left border of the patient image and the right end is on the right border of the patient image. If the line is in the correct position, then we will get the LAT diameter value in the LAT box. In this example, we get the LAT diameter as 31.31 cm.
- Press the **AP** button, a line will appear on the patient image in the AP direction. We have to shift the line and make sure the top edge is on the top border of the patient image and the bottom edge is on the bottom edge of the patient image. If the line is in the correct position, then we will get the AP diameter value in the **AP** box. In this example, we get the diameter of the AP as 22.23 cm.
- Press the **Calculate** button.

We will get the D_{eff} value in the **Diameter (cm)** box. In this example we get the D_{eff} value of 26.38 cm.

Both lines can be removed by pressing the **Clear** button. If desired, the process can be repeated again.

The weakness of this method is that it is very subjective because it depends on the user subjectivity to determine the patient limit. In this case, the appropriate window selection can

be used. In certain cases, the **Fat window** may be better, and in other certain cases the **Bone window** may be better.



Figure 73. Display for calculating D_{eff} manually from a patient image.

b. Automated calculation

For automatic calculation, the user just presses one button (Auto), and IndoseCT will calculate the D_{eff} value automatically. There are 3 options for automatic, namely: Area, Center, and Max (Figure 74). Area is the D_{eff} value calculated from the area of the patient image that has been segmented automatically. Center is the D_{eff} value will be calculated based on the AP and LAT values when passing the patient's midpoint. Max means that the D_{eff} value is calculated based on the maximum AP and maximum LAT values (Anam et al. Atom Indonesia. 2017; 43(1): 55-60).

Calculation the value of D_{eff} with the **Area** method is as follows.

- Select Auto.
- Select Area.
- Press the **Calculate** Button.
- Wait a few moments, then the image segmentation results will be visible and the calculation results for the D_{eff} value will be displayed in the **Diameter (cm)** box (Figure 75). In this case, we get an effective diameter value of 28.05 cm.
- This process wills likely take less than a second.

G IndoseCT v20b.22 □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □		- 🗇 × Phantom: BODY v
	ID Name Age 39 Scanner Sex M TDIvol Diameter SSDE Organ Au Based on: Effective Diameter (Deff) Source: Get from Image Method: Auto Calculate Diameter 0 cm Show Graph	ate 28/01/2019 The second sec
0 100 200 300 400	Previous 500	Next
Sort Images < 24/48 > 24 Go to slice		

Figure 74. D_{eff} calculation with automatic calculation. There are three choices of methods, namely **Area**, **Center**, and **Max**.

In certain cases, this segmentation may not be accurate. In this case, the user can change the threshold value in HU. By default, **IndoseCT** uses a value of -300 HU. In some cases, this threshold results in a fairly good patient segmentation. However, some materials with HU around -200 HU may be included in the segmentation. To overcome this, the threshold value can be increased to -200 HU or other. Under certain conditions, the HU value may need to be lowered, for example to -500 HU.



Figure 75. The display for calculating the D_{eff} automatically from a patient image using the **Area** option. The D_{eff} value obtained is 28.05 cm.

Calculating D_{eff} with the **Center** method is as follows (see **Figure 76**).

- Select Auto.
- Select Center.
- Press the **Calculate** Button.

Wait a few moments, then the image segmentation results will be visible and the calculation results for the D_{eff} value will be displayed in the **Diameter (cm)** box (**Figure 76**). In addition to the **effective diameter** value, we will get values of **LAT** and **AP** diameters. In this example, we get a D_{eff} value of 27.10 cm. This value is generally a little smaller than using the **Area** method.



Figure 76. The display for calculating the D_{eff} automatically from a patient image using the **Center** option. The D_{eff} value obtained is 27.10 cm.

Calculating the D_{eff} value with the **Max** method is as follows (Figure 76):

- Select Auto.
- Select Max.
- Press the **Calculate** button.



Figure 77. The display for calculating the D_{eff} automatically from a patient image using the **Max** option. The D_{eff} value obtained is 27.42 cm.

Wait a few moments, then the image segmentation results will be visible and the calculated D_{eff} value will be displayed in the **Diameter (cm)** box (**Figure 77**). In addition to the D_{eff} value, we also get the **LAT** and **AP** values. In this example, we get a D_{eff} value of 27.42 cm. This value is closer to the **Area** method, compared to using the **Center** method.

The use of effective diameter, reported in several journals, is less accurate for patient dose estimation especially in the chest area. This is because in the chest, there are lungs whose composition is mainly air. The HU value of the lungs is about -600 HU. Thus, the water-equivalent chest diameter should be smaller than the effective diameter. This causes the dose received by the patient to be high.

In order to obtain a more accurate patient diameter and dose estimate, this effective diameter needs to be corrected. The easiest way is to correct the pixels that represent the **Lung**, with a correction factor of 0.3. However, apart from the lungs, there are also bones in the chest with HU values above +1000. **IndoseCT** also provides bone correction options. The **bone** correction factor uses a value of 1.8.

To detect lung, by default **IndoseCT** uses a threshold value of -250 HU, and to detect bone it uses a threshold value of +250 HU. These values need to be tested for accuracy with continuing research. The user can change these threshold values freely.

This lung and bone correction will only be active when we use the **Center** or **Max** options. An example of the **Center** method which uses **Lung** correction only, is shown in **Figure 78**, and one which uses **Lung** and **Bone** corrections is shown in **Figure 79**.

G indoseCT v20b.22	– 0 ×
🖿 🖿 💏 🏟 🖥 🔀 Windowing: Soft Tissue 🗸 360 60	Phantom: BODY ~
Image:	Exam Date 28/01/2019 Institution Scanner SIEMENS-Sensation 64 Protocol BDOMEN Irgan Analyze Options Threshold (HU) _200 © Min. pixel area 500 © Min. pi
Loading Images	

Figure 78. An example of the **Center** method using **Lung** correction. In the image there are 4 color gradations. Outside the patient has a pixel value of 0, lungs have a pixel value of 1, soft tissue has a pixel value of 2, and bones have a pixel value of 3. The D_{eff} value is 21.19 cm.

😋 IndoseCT v20b.22		– 0 ×
🖾 🖿 M 🏟 K Soft Tissue v 360 60		Phantom: BODY ~
100 2	(Deff)	m Date 28/01/2019 itution inner SIEMENS-Sensation 64 ADDOMEN Analyze Options Threshold (HU) -300 © Min. pixel area 500 © Method ○ Area AP 24.69 cm ○ Max Correction Correctio
Previous		Next
A 0 100 200 300 400 500		
Sort Images < 24/48 > 24 Go to slice		

Loading Images

Figure 79. Example of the **Center** method using **Lung** and **Bone** corrections. The D_{eff} value is 22.47 cm.

The effective diameter value obtained is 21.19 cm, which is much smaller than without the **Lung** correction, which is 27.10 cm. The use of **Lung** and **Bone** corrections resulted in an effective diameter value of 22.47 cm, which is slightly larger than using the **Lung** correction alone. The use of **Lung** and **Bone** corrections result in a more accurate patient dose estimate.

c. Z-axis auto calculation

The D_{eff} value is different for each slice. **IndoseCT** gives the user the option to calculate D_{eff} in Z-axis (for all slices). However, the calculation may take a considerable time.

To calculate the D_{eff} of multiple slices, the steps are as follows:

- Select **Diameter** tab.
- Select Effective Diameter (cm) option.
- Select **Get from Image** option. For this, the image must have been opened.
- Select the Auto (Z-axis) option (See Figure 80).
- In this case there are three options: **Slice Step**, **Slice Number**, and **Regional**. The explanation of each will be discussed in the following discussion.

As with the previous **Auto** method, the **Auto** (Z-axis), the **Area**, **Center**, and **Max** options are still applicable. In addition, if the user uses the **Center** or **Max** options, the **Lung** and/or **Bone** corrections are also still valid.



Figure 80. Auto (Z-axis) calculation preview for D_{eff.} In this case, there are three choices, namely **Slice Step**, **Slice Number**, and **Regional**.

The purpose of **Slice Step** is to provide a choice of steps to use. If we select **Slice Step** and fill it with a value of 1, the D_{eff} calculation is done on the first slice, the next slice is the 2nd slice, and the next slice again is the 3rd slice and so on. Because the **Slice Step** we are using is 1, the D_{eff} will be calculated from all slices (**Figure 81**). After we fill the **Slice Step**, we can immediately click **Calculate**.

	1 4 \$ 2 5	Windowin	g: Soft Tissue 🔻	360 60							Phantor	n: BODY
	G Graph of Slice - E	ffective Diameter			-		ID Name Age Sex	39 M		Exam Date Institution Scanner Protocol	28/01/2019 SIEMENS-Sensat ABDOMEN	ion 64
			Slice - Effe	ective Diameter			CTDI	vol Diameter	SSDE Or	gan Anal	yze	
	20	·····		٨	******	******	Bass Effe Sou Get	ed on: ective Diameter (D rce: : from Image	eff) ~	Options Thresho Min. pix Metho	old (HU) -300 🌻 el area 500 🕏 d	
	Deff (cm)						Aut	o (Z-axis) alculate Diamet	v er 27.37 cm	 Are Ce Ma 	ea nter IX	
								Show Graph		 Z-axis C Slice Slice Regi 	options Step Number onal	
	0 Options	10	20	30 slice	40 Save Plo	5 t Close				1]	
					_		Pre	vious			[Next
-100	0	100	200	300	400	500	600					

Figure 81. The D_{eff} calculation display for the Slice Step option equal to 1. The Diameter (cm) box will be filled with the average D_{eff} .

To speed up the calculation of D_{eff} , users can also use **Slice Step** 2 or other values according to the number of slices of patients. The user is advised to use a number of iterations of around 10. For example, if the number of slices is 500, the first to be counted is slice 1, then slice 11, then slice 21, and so on. By filling the **Slice Step** equal to 10, the computation time becomes around 10 times faster than using a **Slice Step** of 1.

In addition to the **Slice Step** option, there is another option to speed up computing, namely using the **Slice Number** option. If we fill the **Slice Number** with 1, then from all existing slices, only 1 slice will be counted in the middle. If we use **Slice Number** 3, D_{eff} will be calculated from 3 slices, namely the middle slice, the third slice from the left and the third slice from the right. **Figure 82** shows the choice of **Slice Number 9**, where the D_{eff} calculation is performed on 9 slices that are evenly distributed along the longitudinal axis. The slice number used to calculate the D_{eff} is shown on the graph.

Another **Auto (Z-axis)** option for calculating D_{eff} is **Regional**. In this case, user can specify a certain slice number to start (e.g. slice number 15) and specify a slice number to end (e.g. slice number 30). This **Regional** option is very useful for calculating the average D_{eff} and SSDE averages for certain organs only.



Loading Images

Figure 82. D_{eff} calculation for the **Slice Number** option. In this case, the **Slice Number** is filled with the number 9. The graph shows 9 points of D_{eff} value. The **Diameter (cm)** box will be filled with the average D_{eff} of these 9 values.

It has been reported by several studies that the organ dose is less accurate if estimated from the whole slice, but it is more accurate if it is estimated only on certain slices where the organ is located (Anam et al. J Biomed Phys Eng. 2021). An example of a Regional option for calculating D_{eff} is shown in **Figure 83**.



Figure 83. D_{eff} calculation display for **Regional** options. In this case, user can specify a specific slice number to start with and a slice number to end. The **Diameter (cm)** box will be filled with the average D_{eff} of the several slices used.

Note: There is currently no protocol for the number of slices to be used for SSDE calculation. The AAPM stated that the calculation of the slice in the middle (along the longitudinal axis) is acceptable. Using **Auto (Z-axis)** option will give a more accurate result, but it will take a longer time.

In the **Auto (Z-axis)** option, if we want the graph to be displayed, **Show Graph** must be selected. As before, this graph can be processed, saved, or exported to Microsoft Excel (**Figure 84**).



Figure 84. The graph of the result of the D_{eff} calculation using **Auto (Z-axis)** can be saved in a certain format.

VII. CALCULATION OF THE WATER EQUIVALENT DIAMETER (D_w)

The water-equivalent diameter (D_w) is a metric that takes into account the size and composition of the patient. Therefore, for patient dose estimation, D_w is more accurate than D_{eff} . **IndoseCT** can be used to calculate D_w . In **IndoseCT**, the user can either directly enter the D_w value because it is already known(**Input manually**), or the D_w value is calculated from the patient image (**Get from Image**) (**Figure 85**).

Steps for D_w calculation are as follows:

- Select the **Diameter** tab.
- Select the Water Equivalent Diameter (D_w) option.
- Choose whether to enter the D_w value manually (Input Manually) or calculate it from the image (Get from Image).



Loading Images

Figure 85. Calculation of D_w there are two choices of methods, namely: Manually and Get From Image. To use the Get From Image option, the user must have opened the patient's CT image in DICOM format.

7.1 Input D_w value manually

The D_w value can be entered manually, by performing the following steps (see Figure 86):

- Select the **Input Manually** option.
- Enter the value of D_w in the **Diameter** field.

G IndoseCT v20b.22		– o ×
🔄 🗈 🥦 👯 🔅 🚺 Windowing: Soft Tissue 🗸 360 60		Phantom: BODY ~
	D Name JB Name Sex M Sex M CTDIvol Diameter SDE CTDIvol Diameter SDE Source: Input Manually Calculate Diameter 25 Cron Source:	019 S-Sensation 64 EN Next
Sort Images < 24/48 > 24 Go to slice		

Figure 86. Display to input the D_w values manually.

7.2 D_w calculation from image

If the user does not have a D_w value, then the D_w value must be calculated from the patient image. In this case there are three options, namely **Manual**, **Auto**, and **Auto** (**Z**-axis) (See **Figure 87**).

C₀ IndoseCT v20b Image: IndoseCT v20b Image: Im	− □ × Phantom: BODY ✓ ID Exam Date 28/01/2019 Name Institution Scanner Age Scanner SIEMENS-Sensation 64 Sex M ✓ Protocol ABDOMEN CIDIvol Diameter SOE Options Options Options
	Water Equivalent Diameter (DW) Source: Get from Image Min. pixel area Source: Min. pixel area Source: Min. pixel area Source: Min. pixel area Cast from Image Auto (Z-axis) Auto (Z-axis) Manual Show Graph Z-axis Options Slice Step Slice Number Regional
500 -100 100 200 300 400 500 600 Sort Images < 24/48 > 24 Go to slice Loading Images	Previous

Figure 87. Calculation of D_w from the image. There are three choices of methods, namely: Manual, Auto, and Auto (Z-axis). To use the Get From Image option, the user must have opened a CT image in DICOM format.

a. D_w calculation manually

To perform manual D_w calculations from images, there are two tools, namely **Ellipse** and **Polygon** (Figure 88). The first D_w calculation, used in AAPM 220, is by using an **Ellipse** tool.



Figure 88. Calculation of D_w from the image manually can use two tools, namely: **Ellipse** and **Polygon**.

D_w calculation with Ellipse tool

To perform D_w calculations from patient images using the **Ellipse** tool, the following steps are taken:

- Select the **Diameter** tab.
- Select the Water Equivalent Diameter (D_w) option.
- Select the **Get from Image** option.
- Select the **Manual** option.
- Then press the **Ellipse** button.
- After that, an ellipse will appear (Figure 89). The user must enlarge, shrink, rotate or move the ellipse so that the ellipse covers the entire patient image and only a few objects outside the patient enter the ellipse.
- In **Figure 89**, the ellipse contains a small circle and a small diamond. The small circle is for rotating the ellipse and the small diamond is for enlarging and shrinking the ellipse, including changing the scale to a circle or an ellipse. To move the ellipse place the cursor inside the ellipse and drag the ellipse to the proper position. If the cursor is outside the ellipse, the shifted object is not the ellipse, but the patient image.
- Even before pressing the **Calculate** button, the D_w value will already be seen in the **Diameter (cm)** box (**Figure 90**).

G IndeseCT v20b.22 Image:	- Ø Phantom: BODY	×
	D Name 39 Scanner Sex M CTDIvol Diameter SSDE Organ Analyze Based on: Water Equivalent Diameter (Dw) v Source: Get from Image Method: Claculate Diameter 18.48 cm Show Graph	
A -100 0 100 200 300 400 500 600	Previous	j
Loading Images < 24/48 > 24 Go to slice		

Figure 89. Calculating D_w values from patient images using the **Ellipse** tool. The small circle on the ellipse is used for rotating the ellipse and the small diamond is used for zooming-in, zooming-out, and changing the scale of the ellipse.

IndoseCT v20b.22	- D Phantom: BODY
pixel: (472, 42) value: -999.000	D Exam Date [28/01/2019] Name J Institution SiteMENS-Sensation 64 Sex M Protocol SiteMENS-Sensation 64 ABDOMEN CTDIvol Diameter SSDE Organ Analyze Based on: Water Equivalent Diameter (Dw) Source: Get from Image Method: Manual Claculate Diameter 24.65 cm
900 g -100 0 100 200 300 400 500 65	Previous

Sort Images < 24/48 > 24 Go to slice Loading Images

Figure 90. The ellipse has covered the entire image and only a few objects outside the ellipse are included within it. The D_w value can be seen in the **Diameter (cm)** box.

D_w calculation with Polygon tool

To perform D_w calculations from patient images using the **Polygon** tool, several steps are taken as follows:

• Select the **Diameter** tab.

- Select the Water Equivalent Diameter (D_w) option.
- Select the **Get from Image** option.
- Select the **Manual** option.
- Then press the **Polygon** button (Figure 91).
- Move the cursor to the patient's image, after which the cursor will change to a cross shape.
- Click on the edge of the patient image, then move to another point (see **Figure 91**) and so on, until all the patient boundaries have been covered with the red plus signs.
- To finish, double-click near the first plus sign. The points we created will then be connected as shown in **Figure 92**.
- Even before pressing the **Calculate** button, the value of D_w will already be seen in the **Diameter (cm)** box (**Figure 92**).
- If there are some points that are not right on the edge of the patient image, these points can be shifted (See **Figure 93**). Points can also be added, by clicking between two existing points.



Figure 91. Display of the segmentation process for calculating D_w values from patient images using the **Polygon** tool.
IndeseCT v20b.22 Image: Im		- D 2 Phantom: BODY
pixel: (274, 116) value: 55.0000	ID Exam Date 28/(Name 39 Scanner SIE Sex M V Protocol ABD CTDIvol Diameter SDE Organ Analyze Based on: Water Equivalent Diameter (Dw) V Source: Get from Image V Method: Manual V Calculate Diamete 24.67 cm Show Graph	1/2019 IENS-Sensation 64 DMEN
	Previous	Next

Figure 92. Display of the final result of the segmentation process using the **Polygon** tool to calculate the D_w value of the patient image.

DD"#\$	Windowing: Soft Tis	sue ~ 360 60					Phanto	m: BODY
	pixel:	(422, 302) value: 51.0000	ID Nome Age Sex CTDiv Base Wat Sour Get Man	39 M Ol Diameter d on: er Equivalent Dian ce: from Image od: ual culate Diameter Diameter	SSDE Org	Exam Date Institution Scanner Protocol an Anal Options Pol Ell CC	28/01/2019 SIEMENS-Sensa ABDOMEN yze ygon ipse ear	tion 64
								March

Figure 93. Example of two points that are shifted outward from the patient boundary and shifted inward from the patient boundary. For an accurate calculation of D_w , the points should be kept within the patient boundary.

b. D_w calculation automatically (Auto)

To perform D_w calculations automatically from patient images, the following steps are taken (see **Figure 94**):

- Select the **Diameter** tab.
- Select the Water Equivalent Diameter (D_w) option.

- Select the **Get from Image** option.
- Select the **Auto** option.
- Next, there are four fields in IndoseCT: Truncated images, Largest object only, Without ROI, and Remove table (Figure 94). These four fields can be left blank (default), or one or two can be filled in.
- Press the **Calculate** button.
- If the calculation is completed, it is indicated by the appearance of the D_w value in the **Diameter (cm)** box and the patient image will be contoured with a red line.



Figure 94. D_w calculation is carried out automatically. There are four fields: **Truncated images**, **Largest object only, Without ROI**, and **Remove table**.

For the default condition, the user does not need to select the four options available. **IndoseCT** will segment the patient objects, even though there may be several objects (such as chest, one hand or both hands), and the patient table is removed automatically (Anam et al. J Appl Clin Med Phys. 2021;1-11). In this default condition, a new image segmentation method is used. The segmentation process and the D_w calculation is performed by pressing the **Calculate** button. The results are shown in **Figure 95**, which shows that the patient image is well segmented, and the D_w value is displayed in the **Diameter (cm)** box.

	- 0
🔟 🔲 🎁 🙀 🚺 🗾 Windowing: Soft Tissue 🗸 360 60	Phantom: BODY
	ID Exam Date 28/01/2019 Name Institution
	Age 39 Scanner SIEMENS-Sensation 64 Sex M V Protocol ABDOMEN
	CTDIvol Diameter SSDE Organ Analyze
100 -	Based on: Options
	Water Equivalent Diameter (Dw) Threshold (HU) 300 ‡ Source: Min. pixel area 500 ‡ Get from Image Truncated image
200 -	Auto
300	Calculate Diameter 24.21 cm
400 -	
500 -100 0 100 200 300 400 500 600	Previous
Sort Images < 24/48 > 24 Go to slice	

Figure 95. Display of automatic calculation of D_w values from patient images by default (without selecting any of the four options: **Truncated image**, **Largest object only**, **Without ROI**, and **Remove table**).

Selecting the default condition and selecting the **Largest object only** will give the same result if the patient object in the image consists of only one part (as in **Figure 95**). But for more than one object, for example the stomach and hands, the result will be different as in **Figures 96** and **97**.



Figure 96. Segmentation using the **Largest object only** option causes only the largest object to be segmented (the patient's stomach), while the patient's hand (which is separated from the abdomen) is not segmented.

G IndoseCT v20b.22 □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □	- D X Phantom: BODY V
0 - 100 - 200	ID Exam Date 17/11/2014 Name Institution Age 54 Scanner F Protocol CHEST CTDIvol Diameter SSDE Organ Analyze Based on: Options Water Equivalent Diameter (DW) ✓ Source: Get from Image Min. pixel area 500 ♀ Method: Largest object only Auto ✓ Show Graph Ferrore table
	500 600
Sort Images < 323/406 > 1 Go to slice	

Loading Images

Figure 96. Segmentation without **the Largest object only** option (default) results in segmentation of body parts, even though in the image there is more than one part (such as the patient's abdomen and hands).

In **Figure 96**, segmentation using the **Largest object only** option results in only the largest object being segmented (the patient's abdomen), while the patient's hand (which is separated from the abdomen) is not segmented (Anam et al. J Appl Clin Med Phys. 2016;17(4):320-333). In **Figure 97**, the segmentation does not use the **Largest object only** option (default) and all sections are well segmented (both the patient's abdomen and hands). Thus, by not selecting the **Largest object only**, a higher D_w value will be obtained(In the old version of **IndoseCT**, the default condition is to use the **Largest object only**. But in this new version, the default condition is not to use the **Largest object only**).

The **Truncated image** option is used to compensate for a cropped image. For an uncropped image as shown in **Figure 95**, the **Truncated image** option will have no effect. However, for a cropped image as shown in **Figure 98**, using the **Truncated image** option results in a more accurate D_w value. In this case, **IndoseCT** will calculate the truncation percentage (TP), then apply a truncation correction factor based on the TP value. Detailed information can be found in Anam et al. Radiat Prot Dosim. 2017; 175(3): 313-320.

Figure 98 shows the calculation of D_w on a cropped image but does not use the **Truncated image** option. In this case, the D_w value is 19.92 cm. **Figure 99** shows the same image calculated using the **Truncated image** optionwhich gives a value of D_w is 20.32 cm. For a truncated image, using a **Truncated Image** produces a slightly larger D_w value.

G IndoseCT v20b.22		− □ × Phantom: BODY ∨
	ID Exam Date Name Institution Age 74 Scaner Scaner Protocol CTDIvol Diameter SSDE Organ Anal Based on: Water Equivalent Diameter (Dw) Thresh Source: Image Image Image Method: Image Image Image Auto Image Image Image Show Graph Show Graph Rem Image Previous Freevious Image Image	DB/04/2016 SIEMENS-Emotion 6 CHEST hyze cold (HU) -300 C cel area 500 C cated image rest object only how table
Sort Images < 9/63 > 1 Go to slice		

Figure 98. Calculation of D_w on a cropped image, but not using the **Truncated image** option. The value of D_w is 19.92 cm.

IndoseCT v20b.22 Image: Diagram Picture Image: Diagram Picture <thimage: diagrampicture<="" th=""> <thi< th=""><th></th><th>− ⊡ × Phantom: BODY ∨</th></thi<></thimage:>		− ⊡ × Phantom: BODY ∨
	D Examined and the second and the se	Im Date 08/04/2016 Titution SIEMENS-Emotion 6 tocol CHEST Analyze Options Threshold (HU) -300 Min. pixel area 500 C Truncated image Largest object only Without ROI Remove table
-100 0 100 200 300 400 500 600		

Figure 99. Calculation of D_w on a cropped image and using the Truncated image option. The value of D_w is 20.32 cm.

Alternatively, the D_w value can be calculated from the entire image although it is less accurate. With this method, the image area exceeds the patient area. However, the linear attenuation of the air is close to 0 and the HU value is around -1000, so it does not have much effect on the final value of D_w (Anam C, et al. Radiat Prot Dosim. 2019; 185(1): 42–49). If we want to use this option, then in the options section we select **Without ROI** (Figure 100). In our

example the D_w value is 25.55 cm, which is slightly larger than using automatic segmentation, 24.21 cm (Figure 95).



Figure 100. Calculation of D_w Without ROI. The value of D_w is 25.55 cm.

The calculation of D_w **Without ROI** is accurate if there is nothing other than air outside the patient. In reality, however, there are sometimes objects with large linear attenuation, for example the patient table. Therefore, the **Without ROI** option can be combined with the **Remove table** as shown in **Figure 101**.



Figure 101. Calculation of D_w using **Without ROI** and **Remove table**. The value of D_w is 24.81 cm.

It can be seen that the patient table was successfully removed and as a consequence the D_w obtained decreased to 24.81 cm. However, this value is still slightly larger than using the two previous options.

Note: By default, image segmentation uses a threshold value of -300 HU (**Figure 102**). This threshold value is not always appropriate. If this threshold value cannot segment the patient properly, the user can change this threshold value to be higher, for example -200 HU, or smaller to -400 HU. In this segmentation, we used a **Min. pixel area** of 500 pixels. If there are objects with an area of less than 500 pixels, then those object will be ignored in the segmentation process. This value can be changed as needed.



Figure 102. For image segmentation automatically using a **Threshold** of -300 HU and **Min. pixel area** of 500 pixels. These values are not always appropriate for segmenting various cases. If these values cannot segment the patient properly, they can be changed.

The special algorithm in **IndoseCT** for segmenting patients has worked well in hundreds of patients with several CT scanners. However, **IndoseCT** also still has limitations, namely in some cases sometimes it cannot segments images properly. If the user judges that the segmentation area is too large or too small, the user can calculate D_w by using manual calculation, i.e. by using the **Polygon** or **Ellipse** tool options.

c. D_w calculation with auto (Z-axis)

The value of D_w is different for each slice. **IndoseCT** gives user the option to calculate D_w along the Z-axis (for all slices) (Anam et al. J Biomed Phys Eng. 2021). However, this calculation may be time-consuming.

These are the steps to calculate D_w from multiple image slices

• Select the **Diameter** tab.

- Select the Water Equivalent Diameter (cm) option.
- Select the **Get from Image** option. For this, the image must be already open.
- Select the Auto (Z-axis) option (See Figure 103).
- Just as with Auto (Z-axis) in D_{eff} calculation, for D_w calculation there are also three options: Slice Step, Slice Number, and Regional.



Figure 103. Auto (Z-axis) calculation display for D_w. There are threeoptions, namely Slice Step, Slice Number, and Regional.

An example of the results of the Auto (Z-axis) calculation graph with Slice Step 1 is shown in Figure 104. The average D_w value is entered in the Diameter (cm) box.

🕸 🕐 🗈 🗹 Windowing: Soft 1	issue ∨ 360 60			Phantom: BODY
G Graph of Slice - Water Equivalent Diame	ter lice - Water Equivalent Diameter	×	ID Name Age 39 Sex M	Exam Date 28/01/2019 Institution
28 26 24 22 20	****	a de la constance de la constan	CTDIvol Diameter SSDE 0 Based on: Water Equivalent Diameter (Dw) Source:	Options Threshold (HU) -300 🗘
20 18 (E 16 (E 16) 14 (C 12)			Get from Image Method: Auto (Z-axis)	Truncated image Largest object only Without ROI
10 8 6 4			Calculate Diameter 26.52 cm ✓ Show Graph	Remove table Z-axis Options Slice Step Slice Step
	20 30 slice	40 50		Regional
Ομισιισ		Jave not Close	Previous	Next

Figure 104. Example of the results of the Auto (Z-axis) calculation graph with Slice Step 1.

Note that D_w calculation with **Auto (Z-axis)** option, and with **Slice Step, 1** option, is the most accurate method of calculating the patient diameter for calculating patient dose, compared to all the previously described methods. This method is recommended for use in the calculation of patient diameter. However, the calculation this may take a longer depending on the number of slices. However, for the number of slices of about 500, the time required is generally not more than 1 minute.

VIII. CALCULATION OF SIZE-SPECIFIC DOSE ESTIMATE (SSDE)

The SSDE value can be calculated after the $CTDI_{vol}$, D_{eff} or D_w values are filled in completely. The SSDE value represents the dose estimate in the individual patient.

To calculate the SSDE value, there are two options, namely calculating the SSDE at one particular slice (**One slice** option) or several slices along the Z-axis (**Z-axis** option).

a. One slice

The steps to calculate SSDE value are as follows (see Figure 105):

- To calculate CTDI_{vol}, first select the **SSDE** tab in the tab rows. If the **SSDE** tab is clicked, the color will change to white.
- Make sure the type of **Phantom** used is **Body** or **Head**.

Note: When calculating the SSDE value, the phantom option must be correctly chosen as **Body** or **Head** phantom, otherwise the SSDE calculation will be wrong. When opening an image, this **Phantom** option will be filled in automatically. However, the user should double-check this.

Next choose the conversion factor to use. For the Body Phantom option, only AAPM 204 is used, while for the Head Phantom, there are two choices, namely AAPM 204 and AAPM 293 (Figure 105). The new and more accurate conversion factor for head is the AAPM 293.

Figure 105. For Body Phantom only conversion factors from AAPM 204 are used, while for Head Phantom there are two choices of conversion factors, namely AAPM 204 and AAPM 293.

• Next, choose the **Protocol** used. For the **Body Phantom** there are 10 choices (**Figure 106**), while for the **Head Phantom** there are three options (**Figure 107**). Actually, the

choice of protocol has no effect on the SSDE value, but it does affect the value of the effective dose. Since the effective dose is calculated in conjunction with the SSDE calculation, the choice of this protocol needs to be done.

- Choose **One slice** option.
- Press the **Calculate** button.





IndoseCT v20b						_	- 0 >
🖿 🖿 🖻 🖫 🏶 🔍 🖬 🗹 Windowing	: Soft Tissue ~ 360 60					P	hantom: HEAD
		TD			Evam Date	20/02/201	6
		Name	, <u> </u>		Institution	20/03/201	0
0 -		Age	67		Scanner	STEMENS-E	motion 6
		Sex	M		× Protocol	HEAD	
100	E	CTDI Base	Ivol Dia r d on:	meter SSDE	Organ Ana	lyze	
		AAP	M 204				~
		Provi I	- Lo -				
		Proto					
		Head	d				~
200 -		Head	d & Neck				
		Neck	ĸ				
		СП	DI _{vol} (mGy)	0.00	DLP (mGy-c	m)	0.00
	Alter and a second	Del	ff (cm)	0.00			0.00
57N			ii (ciii)	0.00	DEr c (moy c		0.00
00		Cor	nv Factor	0.00	Effective Do	ise (mSv)	0.00
1 1/2 -		SSI	DE (mGy)	0.00			
			alculate	Show Graph			
	Contraction of the		-				
00 -			Save				
100 3	Contract of the						
1 Colores	and and the second s						
500							
		Pre	evious				
-100 0 100	200 300 400	500 600					
Cart Images / 6/37 > 1 Co to clico							
Conclandges Color / I Go to sice							
ding Images							

Figure 107. There are three options for selecting the **Protocol** for the **Head Phantom**.

After a moment, the values of $CTDI_{vol}$, D_{eff} (cm) or D_w (cm), Conv Factor, SSDE, DLP, DLPc (corrected DLP) will be displayed. The Conv Factor value is also displayed visually with a graph.

The head phantom conversion factor based on **AAPM 204** is shown in **Figure 108**, and based on **AAPM 293** is shown in **Figure 109**. The conversion factor based on **AAPM 293** is slightly smaller than **AAPM 204**.



Figure 108. Display of CTDIvol, Deff (cm) or D_w (cm) values, Conv Factor, SSDE, DLP, and DLPc (corrected DLP). Conv Factor values can also be displayed graphically. In this case, the AAPM 204 conversion factor is used.



Figure 109. Display of $CTDI_{vol}$, D_{eff} (cm) or D_w (cm) values, Conv Factor, SSDE, DLP, and DLPc (corrected DLP). In this case, the AAPM 293 conversion factor is used.

b. Z-axis

Loading Images

As with $CTDI_{vol}$ and Diameter (D_{eff} or D_w) values, SSDE values can be calculated directly along the Z-axis. To obtain the SSDE value along the Z-axis, it must be ensured that the $CTDI_{vol}$ and Diameter (Deff or Dw) values have also been calculated along the Z-axis. To calculate the SSDE value along the Z-axis, it is done by selecting the **Z-axis** (Figure 110). Next, three options appear as before, namely Slice Step, Slice Number, and Regional. The explanation of these three options has been discussed in the section on calculating $CTDI_{vol}$ and Diameter (D_{eff} or D_w). A graph of SSDE values along the Z-axis can be displayed by selecting Show SSDE Graph (Figure 111).

IndoseCT v20b Im Im	- D X Phantom: BODY ~
D - Name Age Sex	Exam Date 28/01/2019 Institution Scanner 39 Scanner M V Protocol ABDOMEN
100 Based Protoco Chest	on Dameter 3500. Organ Anaryze
200 - Z-2005 Z-2005 Avg.	ce Deff (cm) 27.37 DLP _c (mGy-cm) 674.79
300 Conv Avg.	Factor 1.36 Effective Dose (mSv) 13.40 SSDE (mGy) 13.79 -Z-axis Options culate Show Graph Silce Step O Silce tumber Silce tumber
	ave O Regional
500 -100 0 100 200 300 400 500 600	ous

Figure 110. **Z**-axis mode for calculating SSDE along the Z-axis. There are three further choices, namely **Slice Step**, **Slice Number**, and **Regional**.

G IndoseCT v20b	黎 @ II I Windowing: Soft Tissue ~ 360 60	-
0 -	G Graph of Slice - SSDE — — — X	ID Exam Date Z8/01/2019 Name Institution Institution Age 39 Scanner SIEMENS-Sensation 64 Sex M Protocol ABDOMEN
100 -	Skce - SSDE	CTDIvol Diameter SSDE Organ Analyze Based on: AAPM 204
200 -		Chest v Z-axis v Avg. CTDL _{vit} (m6y) 10.16 DLP (m6y-cm) 497.63
300 -		Avg. Deff (cm) 27.37 DLP _c (m6y-cm) 674.79 Conv Factor 1.36 Effective Dose (mSv) 13.40 Avg. SSDE (m6y) 13.79 Z-axis Options © Slice Step
400 -	2 0 0 10 20 30 40 5 slice	Save Stice Number Regional
-100	Options Save Plot Close	Previous
Sort Images	< 1/48 > 1 Go to slice	

Figure 111. A graph of SSDE along Z-axis for Slice Step 1.

IX. CALCULATION OF EFFECTIVE DOSE

Once the DLP value is obtained it is possible to calculate the patient's effective dose (mSv units). This effective dose reflects the patient's risk and the probability of developing cancer in the future. In the **IndoseCT** software, the effective dose is calculated on an individual patient basis, namely by **DLPc** (corrected DLP), not just DLP-based. To calculate the effective dose, the user must determine the type of protocol used. For **Head Phantom** there are 3 protocol options, namely: **Head**, **Head & Neck**, and **Neck**. For the **Body Phantom**, there are 10 protocol options, namely: **Chest**, **Liver**, **Liver to Kidney**, **Abdomen**, **Adrenal**, **Kidney**, **Chest-Abdomen-Pelvis**, **Abdomen-Pelvis**, **Kidney to Bladder** and **Pelvis**.

The calculation of the effective dose here is coupled with the calculation of the SSDE. Thus, the choices for calculating the SSDE are also used to calculate the patient's effective dose.

The display of the effective dose is shown in **Figure 112**. With the **IndoseCT**, the value of the effective dose is strongly influenced by the size of the patient.





X. SAVE PATIENT DATA

If **CTDIvol**, **Deff (cm)** or **D**_w **(cm)** values, **Conv Factor**, **SSDE**, **DLP**, **DLPc** (corrected DLP) and **Effective Dose** have been calculated, then these values can be saved by pressing the **Save** button (see **Figure 113**). The data is stored in a database, which can be viewed by clicking the **Database** button above the head image (**Figure 113**).



Loading Images

Figure 113. CTDI_{vol}, **D**_{eff} (cm) or **Dw** (cm) values, **Conv Factor**, **SSDE**, **DLP**, **DLPc** (corrected DLP) and **Effective Dose** can be saved by clicking the **Save** button. The data stored in the database is viewed by clicking the button above the image which is given a red box.

The database consists of Patient_ID, Name, Age, Sex, Exam_date, Institution, Manufacturer, Scanner model, Protocol, CTDIvol, Diameter, Diameter type, SSDE, DLP, DLPc, and Effective dose (Figure 114). Patient and institutional names are closed to maintain the confidentiality of patient and hospital data.

Note that **Diameter type** is an explanation of whether the diameter is the effective diameter (D_{eff}) or the water equivalent diameter (D_w) . In the example of **Figure 114**, the **Diameter type** is **Dw**. Also note that the user cannot save both D_{eff} and D_w at the same time. The user can only calculate and save one **Diameter type**. If the user is conducting research and wants to determine the relationship between D_{eff} and D_w , for example, then the data retrieval must be done twice. For example, the D_{eff} data is calculated first, then saved. Then it is repeated again by retrieving D_w data and saving it. The D_{eff} and D_w data are then exported to Microsoft Excel for further analysis.

If the stored data is not visible in the database, the user can press the **Refresh** button. To exit the database, press the **Close** button.

		4 🕸 🖬 🖸	Windowing: Soft Tissue	✓ 360	60									F	hantom:	HEAD
Patie	nts Kecord	1													-	
LAP	id r	atient id	name	906	COV	evam date	institution	manufacturer	model	protocol	CTDIvol	diameter	diameter type	SSDE	DLP	DI Pc
1	iu ș	Judienic_iu	liame	63	F	20141126	institution	TOSHIBA	Aquilion	PELVIS	15	22.8127	Dw	24.0227	484.5	775.93
2				54	м	20141126		TOSHIBA	Aquilion	PELVIS	16	26,7384	Dw	22,1698	572,8	793,68
3				51	F	20141114		TOSHIBA	Aquilion	PELVIS	23,3	28,892	Dw	29,9361	904,04	1161,5
4				53	F	20141114		TOSHIBA	Aquilion	PELVIS	26	29,8231	Dw	32,1969	969,8	1200,9
5				45		20141114		TOSHIBA	Aquilion	PELVIS	15	27,427	Dw	20,1995	589,5	793,84
6				36	F	20141121		TOSHIBA	Aquilion	PELVIS	12,1	26,2466	Dw	17,1293	370,26	524,15
7				49	F	20141121		TOSHIBA	Aquilion	PELVIS	6,1	20,1211	Dw	10,809	229,97	407,49
8				74	м	20141125		TOSHIBA	Aquilion	PELVIS	7,5	22,5182	Dw	12,1454	261	422,66
9				46	м	20141210		TOSHIBA	Aquilion	PELVIS	13,6	25,4925	Dw	19,8084	526,32	766,58
10	D			40	F	20141204		TOSHIBA	Aquilion	PELVIS	11,6	25,8772	Dw	16,6664	460,52	661,65
1	1			59	м	20141209		TOSHIBA	Aquilion	PELVIS	26	29,6074	Dw	32,4811	1190,8	1487,6
12	2			70	F	20141218		TOSHIBA	Aquilion	PELVIS	9,2	24,7759	Dw	13,7136	336,72	501,91
13	3			54	F	20141218		TOSHIBA	Aquilion	PELVIS	13,1	25,7899	Dw	18,8826	482,08	694,88
14	4			45	F	20150109		TOSHIBA	Aquilion	PELVIS	11,6	26,3888	Dw	16,3218	452,4	636,54
15	5			51	F	20150102		TOSHIBA	Aquilion	PELVIS	26	31,056	Dw	30,874	1079	1281,2
10	5			69	F	20150105		TOSHIBA	Aquilion	PELVIS	16,5	27,1036	Dw	22,5071	805,2	1098,3
														_		
ete :	Selected R	.ow(s)												R	efresh	Clo

Figure 114. The database consists of Patient_ID, Name, Age, Sex, Exam Date, Institution, Manufacturer, Scanner model, Protocol, CTDIvol, Diameter, Diameter type, SSDE, DLP, DLPc, and Effective dose.

If there is erroneous or double data (because this database does not use a primary key, so that double data is possible), the user can delete the desired data by clicking on the data to be deleted, then pressing the **Delete Selected Rows** button (**Figure 115**).

	5	~Q ÷				Software In	doseCT 20b_edit1	[Compatibility M	ode] - We	ord			Sign in		-		
File	2	Home Insert I	Design Layout References Maili	ngs	Review	View H	Help 📿 Tell	me what you war	nt to do							A (Shar
<u></u>	×	Cut													P Fin	d v	
G. F	atients	Record													-		>
Ŧ	Export	to Excel															
	id	patient_id	name	age	sex	exam_date	institution	manufacturer	model	protocol	CTDIvol	diameter	diameter_type	SSDE	DLP	DLPC	e
1	1			63	F	20141126		OSHIBA	Aquilion	PELVIS	15	22,8127	Dw	24,0227	484,5	775,933	1
2	2			54	м	20141126		OSHIBA	Aquilion	PELVIS	16	26,7384	Dw	22,1698	572,8	793,68	1
3	3			51	F	20141114		OSHIBA	Aquilion	PELVIS	23,3	28,892	Dw	29,9361	904,04	1161,52	: 2
4	4			53	F	20141114		OSHIBA	Aquilion	PELVIS	26	29,8231	Dw	32,1969	969,8	1200,94	12
5	5			45		20141114		OSHIBA	Aquilion	PELVIS	15	27,427	Dw	20,1995	589,5	793,84	1
6	6			36	F	20141121		OSHIBA	Aquilion	PELVIS	12,1	26,2466	Dw	17,1293	370,26	524,158	1
7	7			49	F	20141121		OSHIBA	Aquilion	PELVIS	6,1	20,1211	Dw	10,809	229,97	407,498	8
8	8			74	м	20141125		OSHIBA	Aquilion	PELVIS	7,5	22,5182	Dw	12,1454	261	422,66	8
9	9			46	м	20141210		OSHIBA	Aquilion	PELVIS	13,6	25,4925	Dw	19,8084	526,32	766,585	i 1
10	10			40	F	20141204		OSHIBA	Aquilion	PELVIS	11,6	25,8772	Dw	16,6664	460,52	661,657	1
11	11			59	м	20141209		OSHIBA	Aquilion	PELVIS	26	29,6074	Dw	32,4811	1190,8	1487,63	12
12	12			70	F	20141218		OSHIBA	Aquilion	PELVIS	9,2	24,7759	Dw	13,7136	336,72	501,918	
13	13			54	F	20141218		OSHIBA	Aquilion	PELVIS	13,1	25,7899	Dw	18,8826	482,08	694,88	1
14	14			45	F	20150109		OSHIBA	Aquilion	PELVIS	11,6	26,3888	Dw	16,3218	452,4	636,549	1
15	15			51	F	20150102		OSHIBA	Aquilion	PELVIS	26	31,056	Dw	30,874	1079	1281,27	2
16	16			69	F	20150105		OSHIBA	Aquilion	PELVIS	16,5	27,1036	Dw	22,5071	805,2	1098,35	; 2,
<	1														-	-	>
Del	ete Se	lected Row(s)												R	efresh.	Clos	e:
-				Abdo	men-F	elvis, Abdom	nen-Pelvis, Kid	ney to Bladde	r and Pe	lvis (lihat Gamb	ar 26, 27						
2				dan 2	8).												

Figure115. The desired data can be deleted by clicking the data, then pressing the **Delete** Selected Rows button.

Furthermore, the data can be exported to Microsoft Excel for further analysis by pressing the **Export to Ecxel** button in the upper left corner (**Figure 116**). To save an Excel file, the user must specify a folder and also a file name. However, this database can also be processed and displayed using **IndoseCT**, without using Excel. An explanation of data analysis and how to display is shown in **Chapter XII**.

C.,			H & B	Windowing: Soft Ticsup y	260	60									p	hantom:	UEAD	
=			48	Windowing. Sole Habae	500	00	1									iunconn.	TIERO	
G, P	Patier	nts Rea	cord													-		-
Ŧ	Ехро	ort to	Excel															
		id	patient_id	name	age	sex	exam_date	institution	manufacturer	model	protocol	CTDIvol	diameter	diameter_type	SSDE	DLP	DLPc	T
1	1				63	F	20141126		TOSHIBA	Aquilion	PELVIS	15	22,8127	Dw	24,0227	484,5	775,933	1
2	2				54	м	20141126		TOSHIBA	Aquilion	PELVIS	16	26,7384	Dw	22,1698	572,8	793,68	1
3	3				51	F	20141114		TOSHIBA	Aquilion	PELVIS	23,3	28,892	Dw	29,9361	904,04	1161,52	2
4	4				53	F	20141114		TOSHIBA	Aquilion	PELVIS	26	29,8231	Dw	32,1969	969,8	1200,94	2
5	5				45		20141114		TOSHIBA	Aquilion	PELVIS	15	27,427	Dw	20,1995	589,5	793,84	1
6	6				36	F	20141121		TOSHIBA	Aquilion	PELVIS	12,1	26,2466	Dw	17,1293	370,26	524,158	1
7	7				49	F	20141121		TOSHIBA	Aquilion	PELVIS	6,1	20,1211	Dw	10,809	229,97	407,498	٤
8	8				74	м	20141125		TOSHIBA	Aquilion	PELVIS	7,5	22,5182	Dw	12,1454	261	422,66	٤
9	9				46	м	20141210		TOSHIBA	Aquilion	PELVIS	13,6	25,4925	Dw	19,8084	526,32	766,585	1
10	10				40	F	20141204		TOSHIBA	Aquilion	PELVIS	11,6	25,8772	Dw	16,6664	460,52	661,657	1
11	11				59	м	20141209		TOSHIBA	Aquilion	PELVIS	26	29,6074	Dw	32,4811	1190,8	1487,63	2
12	12	2			70	F	20141218		TOSHIBA	Aquilion	PELVIS	9,2	24,7759	Dw	13,7136	336,72	501,918	1
13	13	3			54	F	20141218		TOSHIBA	Aquilion	PELVIS	13,1	25,7899	Dw	18,8826	482,08	694,88	1
14	14				45	F	20150109		TOSHIBA	Aquilion	PELVIS	11,6	26,3888	Dw	16,3218	452,4	636,549	1
15	15	5			51	F	20150102		TOSHIBA	Aquilion	PELVIS	26	31,056	Dw	30,874	1079	1281,27	ż
16	16	5			69	F	20150105		TOSHIBA	Aquilion	PELVIS	16,5	27,1036	Dw	22,5071	805,2	1098,35	2
<																		>
Dele	lete S	Select	ed Row(s)												R	efresh	Clos	e

Figure 116. Data can be exported to Microsoft Excel for further analysis by pressing the **Export to Excel** button in the upper left corner.

XI. CALCULATION OF ORGAN DOSE

In addition to CTDI_{vol} , DLP, D_{eff} or D_w , SSDE and effective dose values, **IndoseCT** can also be used to calculate organ dose. The calculation of the organ dose using **IndoseCT** is also very much determined by the diameter of the patient, either D_{eff} or D_w . Since the dose of the organ is strongly influenced by the patient's diameter and SSDEthe CTDIvol, diameter (D_{eff} or D_w), and SSDE values must be calculated beforehand.

To perform an organ dose calculation, the **Organ** tab must be selected first. There are two choices of methods for calculating organ doses, namely **MC (Monte Carlo) Data** and **Direct Calculation**, as shown in **Figure 117**.

However, it should be understood that the calculation of the dose of this organ still has many limitations and challenges. The accuracy of the organ dose using IndoseCT needs further verification.



Sort Images < 84/216 > 1 Go to slice

Figure 117. There are two choices of methods for calculating organ doses, namely MC (Monte Carlo) Data and Direct Calculation.

10.1 Monte Carlo data

Calculation of organ dose with Monte Carlo data (**MC Data**) refers to the results of a previous publication (Sahbaee et al. Med Phys. 2014; 41(7): 072104). This research was conducted using a Monte Carlo simulation on computational phantoms of varying sizes. From the study, the relationship between SSDE and dose in several organs was determined for several protocols. This calculation method is only applicable to the available protocols. The diameter used in Sahbaee et al (2014) only uses the effective diameter. In **IndoseCT**, this approach is extended to the water equivalent diameter (D_w).

There are 28 sensitive organs estimated in this software, namely: Marrow, Bones, Skin, Brain, Eyes, Larynx-Pharynx, Thyroid, Trachea-Bronchi, Esophagus, Lungs, Thymus, Breasts, Heart,

Liver, Stomach, Spleen, Large Intestine, Adrenals, Pancreas, Small Intestines, Kidneys, Gallbladder, Ovaries, Uterus, Vagina, Bladder, Prostate, and Testes.

The steps for estimating organ dose values are as follows:

- Previously, the CTDIvol, diameter, and SSDE values must be obtained.
- Select the **Organ** tab. The color will change to white.
- Next select the protocol used. There are 13 protocol options. For the head there are three protocols. namely: Head, Head & Neck, and Neck. For the body there are ten protocols, namely: Chest, Liver, Liver to Kidney, Stomach, Adrenal, Kidney, Chest-Stomach-Pelvis, Stomach-Pelvis, Kidney to Bladder and Pelvis (see Figure 118).
- Press the **Calculate** button. The organ dose values will be obtained and displayed on the organ dose boxes (**Figure 118**) and are also displayed visually with a graph (**Figure 119**).



Figure 118. An example of calculating the organ dose in the **Head** protocol using the **MC Data** method.

As with previous graphs, this organ dose graph can be saved or exported to Microsoft Excel. It can be seen from **Figure 118** that the radiation dose for organs located in the head area are large (In this example, the dose in **Brain** is 47.12 mGy) because it is exposed to the primary radiation. For the organs that are outside the head (e.g. **Pancreas**, **Uterus** and so on), the radiation dose is small, because they are only exposed to scattered radiation.



Figure 119. Display graph for organ doses. The x-axis is the initial of the organ name. Organ names can be seen on the main **IndoseCT** form as shown in **Figure 118**.

10.2 Direct calculation

The dose to the organ can also be estimated directly from the patient image (Refer to a previous study on dose distribution between the periphery and center of phantoms of varying size (Anam et al. J Xray Sci Med. 2020; 28: 695-708)). The dose at a point (organ) that is at a certain position is calculated by linear interpolation between the dose in the center (called as the SSDE central, SSDE_c) and the dose at the edge (called as SSDE peripheral, SSDE_p). SSDE_c and SSDE_c are calculated from SSDE_w using k and h correction factors (*k*-factor and *h*-factor).

The organ dose in **IndoseCT** is not only calculated at one point (one pixel), but is calculated from all the pixels in the organ. For this purpose, the user must perform organ segmentation manually. From here, the effective distance map (DEM) will be obtained if using D_{eff} , or the water-equivalent distance map (DWM) if using D_w . Each point on the DEM or DWM is then calculated using linear interpolation to give a dose map (DM). The average DM value from all pixel values within the segmentated area is calculated and this value indicates the estimated value of the organ dose.

The steps for calculating the organ dose using the **Direct Calculation** method are shown in **Figure 120**.

- Previously, the CTDI_{vol}, diameter, and SSDE values must be obtained.
- Select the **Organ** tab. The color will change to white.
- Select the **Direct Calculation** option.
- Click Add contour. When the cursor is hovered over the image, the cursor changes from an arrow to a positive sign.
- Then make a segmentation of the selected organ. To do this, click on the edge of the organ several times so that it covers the organ (Figure 121). At the last point, double-click, then the points will be connected to each other (Figure 122).

• Click the **Calculate** button. The organ dose will be calculated. The mean dose value is shown in the **Mean (mGy)** box and the standard deviation is shown in the **Std**. **Deviation** (mGy) box.



Figure 120. The dose of the organ is calculated directly by selecting the **Direct Calculation** option and followed by making a contour with **Add Contour**.

🔓 IndoseCT v20b.25			– o ×
🖿 🖿 🖷 📽 👰 🗊 🗹 Windowing: Soft	Tissue ~ 360 60		Phantom: BODY V
pixel: (1	81, 366) value: -50.0000	ID Exam Date 2 Name Institution Age 39 Scanner Sex M Protocol CTDIvol Diameter SSDE Organ Analyze Method: Direct calculation Deff (cm) 28.05 Mean (m6y) SSDE_ (m6y) 5.59 Std. Deviation Std. Deviation SSDE_ (m6y) 0.00 Std. Deviation	8/01/2019 EMENS-Sensation 64 BDOMEN a a (mGy) 0.00
400 -		Options Quick Mo Clear Contour	de active Distance Map
500 A) -100 0 100 2	00 <u>300</u> 400 500	Calculate Show Do	se Map togram rr. Factor Graph
Sort Images < 24/48 > 24 Go to slice			
Type here to search	o 🛱 💽 🐂 😭 🚿 🚱	🕼 25°C Hujan ringan \land 🖗 🖩	⊃ /‰ ⊄») 5:49 07/09/2021 ₹ 3

Figure 121. Segmentation of the organ for which the dose will be calculated by clicking several times on the edge of the organ so that it covers the organ.

In **Figure 122**, not only the average organ dose value is shown, but the values of **SSDEw**, **SSDEc** and **SSDEp** are also shown. The symbol **w** refers to weighted, i.e. the average SSDE (weighted)

obtained from $CTDI_{vol} \times f$. The average organ dose calculated on the basis of D_{eff} in this example is 7.16 mGy.

G IndoseCT v20b.25	– 🔿 🗙
🔤 🖿 🖿 🛱 🏟 🕐 💽 🖍 Windowing: Soft Tissue 🗸 360 60	Phantom: BODY ~
	ID Exam Date 28/01/2019 Name Jnstitution Scanner Sex M Protocol ABDOMEN CTDIvol Diameter SSDE Organ Analyze Method: Direct calculation Image: CTDLyoid (mGy) (7.16) Deff (cm) 28.05 Hean (mGy) 7.16 SSDE, (mGy) 5.59 Std. Deviation (mGy) 0.50 SSDE, (mGy) 5.69 SSDE, (mGy) 0.50 SSDE, (mGy) 8.24 Options Quick Mode Clear Contour Show Dise Map Show Histogram Show Korr. Factor Graph
Sort Images < 24/48 > 24 Go to slice Loading Images	
■ P Type here to search O H C = C = 0 G = 0	🌾 25°C Hujan ringan ∧ টু 🗈 /⁄k 40) 5:50 🔖

Figure 122. If the contour or segmentation is completed, the organ dose is calculated by pressing the **Calculate** button. The mean value and standard deviation of the organ dose will be calculated. This calculation is based on $D_{eff.}$ In this example, the mean value of the organ dose is 7.16 mGy.

If we want to display the *h*-factor and *k*-factor, dose map (DM), and histogram dose map (HDM), then these parameters must be selected before pressing **Calculate**, as shown in **Figure 123**. An example of *h*-factor and *k*-factor is shown in **Figure 124**. An example of a dose map (DM) and a histogram of dose map (HDM) for the D_{eff} basis is shown in **Figure 125**. The dose map is radial because it uses a D_{eff} basis.

The organ dose calculated using D_w is shown in **Figure 126**. The average organ dose is 7.77 mGy. This value is slightly larger than using D_{eff} . An example of a dose map (DM) and a histogram of dose map (HDM) based on D_w is shown in **Figure 127**. The dose map is no longer radial when using D_w .

G IndoseCT v20b.25		– o ×
🖿 🖿 📬 🏟 🔍 🖬 🖸 Windowing: Soft	Tissue ~ 360 60	Phantom: BODY V
0-		ID Exam Date 28/01/2019 Name Institution Institution Age 39 Scanner SEMENIS-Sensation 64 Sex M Protocol ABDOMEN CTDIvol Diameter SSDE Organ Analyze Method: Direct calculation V Deff (cm) 28.05 Hosp (mEr) 7.15
200 -		CTDL _{val} (mGy) 5.59 Std. Deviation (mGy) 0.50 SSDE _w (mGy) 7.39 .509
400 -		Options Quick Mode Class Centour Show Effective Distance Map
500 - -100 0 100 20		Calculate Show Vose Map
Sort Images < 24/48 > 24 Go to slice		
Type here to search	o # C 듺 🕿 🕖 🖬	

Figure 123. If the user wants to display the h-factor and k-factor, dose map, and histogram dose map, these parameters must be selected before pressing **Calculate**.



Figure 124. Example view of *h*- and *k*-factors.



Figure 125. An example of a dose map (DM) and an organ histogram of the dose calculated using D_{eff} .

G IndoseCT v20b.25		- <u> </u>
🖿 🖿 🐚 📆 🕸 📿 📑 🚺 Windowing: Soft	Tissue v 360 60	Phantom: BODY ~
0 - 100 - 200 - 300 -		ID Exam Date 28/01/2019 Name Institution Institution Sex M Protocol SEMENS-Sensation 64 ABDOMEN SIEMENS-Sensation 64 ABDOMEN CTDIvol Diameter SSDE Organ Analyze Method:
	200 300 400 500 600	Options Quick Mode Show Water Eq. Distance Map Calculate Show Dose Map Show Histogram Show Corr. Factor Graph
Cont Images < 24/48 > 24 Go to slice		
Type here to search	o # C ≓ 🕿 🕖 #	◇ 23°C Hujan ringan へ Ĝ E

Figure 126. Mean value and standard deviation of organ dose calculated using D_w . In this example, the mean value of the organ dose is 7.77 mGy.



Figure 127. An example of a dose map and an organ dose histogram calculated using D_w.

XII. DATA ANALYSIS

After the patient data is stored, the data can be analyzed to provide useful information or reports to the authorities. The results of the analysis can be displayed in graphs. Reporting and documenting patient doses is very important, although sometimes it has not received much attention from the authorities.

To activate the data analysis section, click the **Analyze** button. The button will change to white, and a form like **Figure 128** will be displayed

In analyzing the data, user can filter data based on: Institution, Manufacturer, Scanner, **Protocol, Sex, Age**, and **Exam Date**. The data entry for each filter is based on the data that has been stored in the database. For example, **Figure 128** shows the protocols of the stored data, namely **PELVIS, HEAD, CHEST, ABDOMEN, CHEST_TO_PELVIS, LEG, Unspecified**, and **All**. **Unspecified** means that the data is not available (filled in) or information has been retrieved from DICOM info, but no information is available in the protocol section.



Figure 128. Display for analyzing patient dosimetry data.

Information on the number of patients appears at the bottom, namely the **Data Count**. In the example of **Figure 129**, no filtering is performed (so the filter is **All**), except for the **Protocol** which uses the **CHEST** option. It is shown that the number of patients is 137.

The user can select the x-axis and y-axis (Figure 129). On the x-axis there are several options, namely Record_ID, CTDI_{vol}, Age, Deff, Dw, SSDE, Effective Dose, DLP, DLPc, and Frequency (see Figure 130). On the y-axis there are also the same options (see Figure 131).

Ga IndoseCT v20b.23	Phantom:	ال HEAD	×
0 - Name Exam Date 10/0 Name Institution Age - Scaner Man Sex Unspecified Protocol CTDIvol Diameter SSDE Organ Analyze	8/2021 ufacturer-Mode	el	
200 - CHEST CHART	s selection Ixis Dw Ixis CTDIvol Genera	~ ~ ate]
300 - Sex P Image: Control of the image: Control of th			
500 -			
600 └ _100 0 100 200 300 400 500 600 700			

Figure 129. After filtering, the x-axis and y-axis must be determined.

G IndoseCT v20b.23		- 0 ×
🖾 🖿 📭 🛱 🕸 🕐 🖥 🗹 Windowing: Soft Tissue 🗸 360 🛛 60		Phantom: HEAD \sim
0 -	ID Exam Date 11 Name Institution Age - Scanner M Sex Unspecified V Protocol CTDIvol Diameter SSDE Organ Analyzi)/08/2021 anufacturer-Model
200 -	Filter Institution All Scanner	Axis selection x-axis Record ID y-axis Record ID CTDIvol Age Deff Dw SDE
300 -	Sex All Image: Constraint of the second sec	Effective Dose DLP DLPc
	8	
Sort Images < 0/0 > 0 Go to slice READY		

Figure 130. There are options on the x-axis, namely Record_ID, CTDI_{vol}, Age, D_{eff}, D_w, SSDE, Effective Dose, DLP, DLPc, and Frequency.

G IndoseCT v20b.23			– 0 ×
🖾 🖿 🐚 🖫 🏟 🕐 💽 🔀 Windowing: Soft Tissue 🗸 36	0 60		Phantom: HEAD ~
0 -	ID Name Age Sex		Exam Date 10/08/2021 Institution Scanner Manufacturer-Model Protocol
100	Filte	vol Diameter SSDE Org: r itution All	Avis selection v axis Record ID v
200 -	Man Scar Prot	ufacturer All nner All ocol All	y-axis CTDIvol ✓ CTDIvol Age Deff Dw
300	Sex Age Exar Date	All -	SSDE Effective Dose DLP DLPc Frequency
400 -		Reset Filter	
500			
-100 0 100 200 300	400 500 600 700		

Figure 131. There are options on the y-axis, namely Record_ID, CTDI_{vol}, Age, D_{eff}, D_w, SSDE, Effective Dose, DLP, DLPc, and Frequency.

Figures 132-139 are examples of some of the graphs of patient data for analysis. Figure 132 is a graph of Record_ID vs CTDI_{vol} using the CHEST protocol. The random data do not form a certain pattern. Figure 133 is a graph of D_w vs CTDI_{vol} (using the same protocol). The value of CTDI_{vol} increases with an increase in the value of D_w , up to a certain D_w value. The increase in the CTDI_{vol} value was due to the use of the tube current modulation (TCM) method.



Figure 132. Example of a graph between Record_ID vs CTDIvol using the CHEST protocol.



Figure 133. An example of a graph between D_w vs CTDI_{vol} using the CHEST protocol.

The graph in **Figure 133** can be adjusted in size. An example of an enlarged graph is shown in **Figure 134**. Trendlines can also be added to the graph. In the example of **Figure 134**, the added **Trendline** is **Linear**. A red straight line appears with the equation and R^2 value of the relationship between the two quantities **D**_w and **CTDIvol.**



Figure 134. D_w vs CTDI_{vol} relationship graph which has been enlarged and given a Linear Trendline.

The **trendline** on the graph can be easily replaced and adjusted according to the data obtained. Figure 135 is an example of a D_w vs CTDIvol relationship graph that has been

enlarged and given a **Trendline** in the form of a **2nd order Polynomial**. Besides **Linear** and **Polynomial** trendlines, there are also **Exponential** and **Logarithmic**. To remove the **Trendlin**, press the **None** option. To display the x-axis and y-axis scaling adjustments by right-clicking on the graph, options will appear as shown in **Figure 135**. Highlight **X-axis**, then fill in the minimum and maximum values. Repeat for the y-axis as needed. Several options exist, for example **Plot Option** and **Export**.



Figure 135. Dw vs CTDI_{vol} relationship graph that has been enlarged and given a Trendline in the form of a **2nd order Polynomial**. The x-axis and y-axis scales have been adjusted.

In addition to displaying the trendline of the data, **lines** can also be added to the graph showing the average value and standard deviation, either on the x-axis or on the y-axis. **Figure 136** is the same graph as **Figure 135**, but with the addition of lines of mean values and standard deviations on the x-axis and y-axis.

IndoseCT v20b.23	Windowing: Soft Tissue V 360 60	- 🗇 > Phantom: HEAD
Options ? X Mean X-data y-data Standard Deviation X-data y-data Trendline None Linear Polynomial Close 500	Graph of Dw - CTDIvol	ID Exam Date 10/08/2021 Institution Institution Sex Unspecified Protocol CTDIvol Diameter SSDE Organ Analyze Filter Analyze Filter Institution All ✓ Y: +3.44 mGy Protocol CHEST ✓ Y: -3.44 mGy Sex All ✓ 40 ave Plot Close Institution Isseet Filter
500 -100 0	100 200 300 400 500	600 700
Sort Images < 0/0 >	0 Go to slice	

Figure 136. An enlarged graph of D_w vs **CTDI**_{vol} with the addition of the mean and standard deviation lines on the x-axis and y-axis.

Figure 137 shows the relationship between D_w vs $CTDI_{vol}$ (same as Figures 134-136 but with additional filters). In addition to the **Protocol** filter, a **Sex** filter is added, which is for women (**Female, F**). The **Data Count** only 137 patients. Thus, we can directly compare the relationship of D_w and $CTDI_{vol}$ between male and female patients.



Figure 137. Graph of D_w vs $CTDI_{vol}$ for data filtered by the CHEST protocol and female gender (F).

Figure 138 is a graph between D_w vs SSDE (using the same protocol) with a Trendline of 2^{nd} Order Polynomial. The SSDE value increases along with the increase in the value of D_w to a certain point, after which it decreases.



Figure 138. An example of a graph of **D**_w vs **SSDE** using the CHEST protocol.

Figure 139 is a graph of Age vs D_w (using the same protocol) using a Linear Trendline. The value of D_w tends to decrease as Age increases.



Figure 139. A graph of Age vs Dw using the CHEST protocol.

Figure 140 is a graph of D_w vs **Frequency** (using the same protocol). Since the y-axis is **Frequency**, this is a histogram of D_w . The user can specify the number of bins to form this histogram. In this example the number of bins is 15. In addition, the user can add a line showing the mean and standard deviation of the histogram.



Figure 140. Example of a graph of D_w vs Frequency.

Figure 141 is a graph of **Age** vs **Frequency** (using the same protocol). The mean age of the patients is about 50 years.



Figure 141. Example of a graph of Age vs Frequency.

The graph of CTDI_{vol} vs Frequency (using the same protocol) is shown in Figure 142. The graph of SSDE vs Frequency is shown in Figure 143. The graph of Effective dose vs Frequency is shown in Figure 144.







Figure 143. Example of a graph of SSDE vs Frequency.

	🕸 🕐 🖪 Ζ Windowing: Soft Tissue → 360 60			Phantom: HEAD
	Graph of Effective Dose - Frequency	ID Name Age -	Exam Da Institutio Scanner	n Manufacturer-Model
	Effective Dose - Frequency	CTDIvol	Diameter SSDE Organ Ar	nalvze
		Filter		Axis selection
	40 Frequency	Institutio	on All ~	x-axis Effective Dose ~
		Manufac	cturer All v	y-axis Frequency V
	30	Scanner	r All V	Bins 15
	neuch	Sex	All	Generate
	<u>ل</u> ۲۵ ک	Age	- 🗘 to - 🗘	
		Exam Da	ate 01/11/2014 🔹 to 19/11/2015 🖨	
		Data Co	Reset Filter	1
			Nuber Filter	1
	₀⊢───└╼─└╼─			
	8 10 12 14 16 18 20 22 24 26 28 30	32 34 36 38		
	Options S	ave Plot Close		
-100	0 100 200 300 400 500	600 700		
ort Images	< 0/0 > 0 Go to slice			

Figure 144. Example of a graph of effective dose vs Frequency.

Figures 140-144 are histograms of $CTDI_{vol}$, SSDE and effective dose values. From these histograms, we can see the distribution of patient doses. This data can be used to determine the diagnostic reference level (DRL).
XV. SUMMARY

IndoseCT 20.b is software for calculating individual patient doses, storing the data in a database, and analyzing the data. With **IndoseCT**, the patient's dose can be controlled in a simple way. With this software, individual patient doses will be known and their comparison to the population can be easily shown.

This software is a tool to monitor radiation doses of all patients undergoing CT scans. If there is a human error in the CT application, then it will be easily monitored. Special strategies can also be designed if the patient dose at the hospital is higher than it should be.

XVI. REFERENCES

- 1. Andersson J, Pavlicek W, Al-Senan R, et al. Estimating patient organ dose with computed tomography: A review of present methodology and required DICOM information. A Joint Report of AAPM Task Group 246 and the European Federation of Organizations for Medical Physics (EFOMP). College Park, MD: American Association of Physicists in Medicine; 2019.
- 2. Ban N, Takahashi F, Sato K, et al. Development of a web-based CT dose calculator: WAZA-ARI. Radiat.Prot Dosimetry. 2011;147:333–337.
- 3. Boone JM, Strauss KJ, Cody DD, et al. Size-specific Dose Estimates (SSDE) in Pediatric and Adult Body CT Examinations: Report of AAPM Task Group 204. College Park, MD: American Association of Physicists in Medicine; 2011.
- 4. Boone JM, Strauss KJ, Hernandez AM, et al. Size-specific dose estimate (SSDE) for head CT: Report of AAPM Task Group 293. College Park, MD: American Association of Physicists in Medicine; 2019.
- 5. Dixon RL, Anderson JA, Bakalyar DM, et al. Comprehensive Methodology for the Evaluation of Radiation Dose in X-ray Computed Tomography: A New Measurement Paradigm Based on a Unified Theory for Axial, Helical, Fan-beam, and Cone-beam Scanning with or without Longitudinal Translation of the Patient Table. College Park, MD: American Association of Physicists in Medicine; 2010.
- 6. Edyvean S. CT equipment and performance issues: radiation protection 162. Radiat Prot Dosimetry. 2013;153(2):190-196.
- 7. Ikuta I, Warden GI, Andriole KP, Khorasani R, Sodickson A. Estimating patient dose from x-ray tube output metrics: Automated measurement of patient size from CT images enables large-scale size-specific dose estimates. Radiology. 2014;270:472-480.
- 8. Leng S, Shiung M, Duan X, Yu L, Zhang Y, McCollough CH. Size-specific dose estimates for chest, abdominal, and pelvic CT: effect of intrapatient variability in water-equivalent diameter. Med Phys. 2015;276(1):184-190.
- 9. McCollough C, Cody D, Edyvean S, et al. The Measurement, Reporting, and Management of Radiation Dose in CT: The Report of AAPM Task Group 23: CT Dosimetry. College Park, MD: American Association of Physicists in Medicine; 2008.
- 10. McCollough C, Bakalyar DM, Bostani M, et al. Use of Water Equivalent Diameter for Calculating Patient Size and Size-Specific Dose Estimates (SSDE) in CT: The Report of AAPM Task Group 220. College Park, MD: American Association of Physicists in Medicine; 2014.
- 11. McMillan K, Bostani M, Cagnon C, Zankl M, Sepahdari AR, McNitt-Gray MF. Size-specific, scanner-independent organ dose estimates in contiguous axial and helical head CT examinations. Med Phys. 2014;41:121909.
- 12. Mihailidis D, Tsapaki V, Tomara P. A simple manual method to estimate waterequivalent diameter for calculating size-specific dose estimate in chest computed tomography. Br J Radiol. 2021;94:20200473.
- 13. Sahbaee P, Segars WP, Samei E. Patient-based estimation of organ dose for a population of 58 adult patients across 13 protocol categories Med. Phys. 2014;47:072104.